

premium
plus™

dental products 2024



innovators, designers, manufacturers,
partners and much, much more...

European Catalogue

INDEX

P.1-2 INTRODUCTION

Who are we
Why choose us
What do we make
Where are we

P.3-8 SMALL EQUIPMENT

Apex Locators
Curing Lights
Curing Light Meter
Low Speed Handpieces

P.9-14 PROPHY

Prophy Angles
Prophy Brushes
Prophy Cups
Prophy Points

P.15-16 Green Line

Disposable Bamboo Cups
Disposable Micro Applicators
Disposable Micro Applicator Tips
Disposable Syringe Tips

P.17 SYRINGE TIPS

Disposable Syringe Tips

P.18 DISPOSABLE

Disposable Bibs
Disposable Face Masks
Disposable Napkin Holders
Disposable Tray Cover Papers

P.19-20 BARRIERS

Barrier Film
Bite Block Sleeves
Chair Sleeves
Curing Light Guide Sleeves
Curing Light Sleeves
Flowable Composite Sleeves
Handpiece & Pen Sleeves
Headrest Covers Non Woven
Headrest Sleeves
Impression Gun Sleeves
Instrument Tray Sleeves
Intraoral Camera Sleeves
Keyboard Sleeves
Mouth Mirror Sleeves
PC Mouse Sleeves
PSP Barrier Envelopes
Suction Tube Sleeves
Syringe Sleeves
T-Style Light Handle Sleeves
Tubing Sleeves
X-Ray Head Sleeves
X-Ray Sensor Sleeves

P.21 APPLICATORS

Brush Applicators
Brush Tips and Handle
Micro Applicators
Micro Applicator Dispensers
Micro Applicator Tips

P.22 EVACUATION

Aspirator Brushes
HVE Suction Tubes
Oral Evacuators
Saliva Ejectors
Suction Tube Adaptors
Surgical Aspirators

P.23-24 IMPRESSION

Alginate and Plaster Spatulas
Autoclavable Impression Trays
Autoclavable Mixing Bowls
Disposable 3-in-1 Trays
Disposable Impression Syringes
Disposable Impression Trays
Disposable Intraoral Syringe Tips
Disposable Mixing Tips

P.25 ENDO INSTRUMENTS

Aluminium Endo Holder
Endo Dispensers and Holders
Endo Measuring Blocks
Endo Rulers
Endo Sponges
Rubber Dam Frame

P.26-31 SUNDRIES

Acrylic Mixing Cups
Aluminium Bur Blocks
Autoclavable Bur/Endo Organisers
Autoclavable Cheek Retractors
Articulating Mounting Plates
Autoclavable Instrument Trays
Autoclavable Mouth Lippers
Autoclavable Napkin Holders
Autoclavable Prophy Rings
Brass Bur Brush
Dental Appliance Boxes
Dental Bib Holder and Clip
Dental Bibs and Cups
Disposable Bite Loops
Disposable Cheek Retractors
Disposable Composite Tips
Disposable Curved Tip Syringes
Disposable Dappen Dishes
Disposable Evacuation Traps
Disposable Eye Shields
Disposable Face Masks
Disposable Fluoride Trays
Disposable Healthy Teeth Bibs
Disposable Healthy Teeth Cups
Disposable Healthy Teeth Masks

P.26-31 SUNDRIES

Disposable Mixing Pads
Disposable Mixing Wells
Disposable Mouth Lippers
Disposable Mouth Mirrors
Disposable Plastic Wedges
Disposable Prophy Rings
Disposable Sterilisation Pouches
Disposable Sticky Applicators
Fracture Detector
Glass Dappen Dishes
Glass Mixing Slab
Jumping Teeth
Interdental Brushes
Magnetic Bur Stands
Mask Extenders
Mouth Props
Needle Station
Patient Tongue Cleaners
Pen/Toothbrush Holder
Plastic Bur Dispensers
Retainer/Mouth Boxes/Denture Baths
Silicone Instrument Mat
Silicone Instrument Rings
Wire Brushes
X-Ray Bite Cushions
X-Ray Film Holders

P.32-33 PLASTIC PRODUCTS

Barrier Film Dispenser
Bib and C-Fold Towel Dispensers
Chain Dispenser
Composite Dispensers and Organisers
Composite Syringe Holders
Cotton Roll Dispenser
Cup Dispensers
Glove and Tissue Dispensers
Grinding Box
Impression Gun Holders
Light Protection Covers
Light Protection Hand Shield
Mask Dispenser
Micro Applicator Dispensers
Pouch Dispensers

P.34-35 ENDO INSTRUMENTS

Barbed Broaches
Gates Drills
Ni-Ti Files
Peeso Reamers
Pluggers and Files
Spiral Fillers
Stainless Steel Files

P.36-39 PIEZO SCALERS

LED Scalers
Scaler Handpieces
Scaler Tips
Scaler Wrenches and Accessories

Who are we ?



At **Premium Plus**, manufacturing dental products is what we do. Only for the past decade we've not just made dental products... we've made better dental products.

Through innovation, research, investment in manufacturing, technical know-how and over 40 years experience in the dental trade, we've learned how to control and improve every step of the manufacturing process... and that's a huge advantage to our customers.


With a global view of the dental industry and knowledge of the different geographical markets, we have learned to understand what you expect and are committed to far exceeding your requirements. Our customers also benefit from our centrally located warehouse in Europe.



Why choose us ?


Whether its under our own **Premium Plus** brand, private labels or OEM projects, our emphasis is on quality for not only the products, but also its presentation. We not only manufacture in-house – We innovate and develop products, research raw material selection, manufacture moulds, construct production lines, create and design artwork and packaging for our partners, warehouse goods and turn our customers' ideas and requests into reality... providing the complete package.

Our continual investment in purpose-built machinery, technical training and market research along with our flexible approach keeps us at the forefront of dental manufacturers.



What do we make?

Over the past decade our range has increased to cover small equipment and the majority of plastic and paper disposable dental products. This gives us, and more importantly you, the ability to offer a huge range that includes LED curing lights, apex locators, handpieces, disposable prophylaxis cups and angles, micro applicators, syringe tips, barrier sleeves, saliva ejectors and suction tubes etc.



disposable oral evacuators

disposable saliva ejectors



LED curing lights



disposable syringe tips



disposable trays for implants



disposable 3-in-1 trays



disposable prophylaxis angles




disposable prophylaxis cups



disposable prophylaxis brushes




disposable micro applicators



disposable PSP barrier envelopes

Where are we?

Premium Plus now supplies to over 90 countries and regions throughout the world and has distribution partners in the USA, Japan, Poland and the United Kingdom. The manufacturing site in China is constantly working on not only maintaining but exceeding standards set by ISO, CE and FDA. We are proud of the fact that we hold **ISO13485: 2016** certification for the **Quality Management System of Medical Devices** and have all of the class IIa products manufactured by us certified with the **CE₀₁₂₃** mark by **TUV SUD**. We know that with our experience, our integrated range of production processes and our ever-expanding range we have the ability to cover all our customers' requirements.



LOW SPEED HANDPIECES AND AIR MOTOR



12 MONTH WARRANTY

A08 DENTAL APEX LOCATOR

8 colours available:



- Bright LCD screen with the trajectory of the file being clearly indicated by different colours
- Foldable design with adjustable view angle, saves space and protects the screen
- The unit is designed with a smooth seamless body and membrane buttons to enhance cold sterilization and help avoid cross contamination
- Advanced multiple frequency technology with automatic calibrating ensures accurate measurement
- Measurement results are not affected by tooth type, size of root canal or the internal environment of root canal
- Different beep tones indicate the position of the file in the root canal
- Screen size: 110 mm x 70 mm Overall size: L122 mm x W118 mm x H25 mm
- Accessories autoclavable up to 134°C
- Weight: 300 gms (including rechargeable lithium battery)

A08 Apex Locator
A08-1 Battery
A08-2 Spare Probe
A08-3 Spare Hook
A08-4 Spare Cord
A08-5 Spare Clip

24 MONTH WARRANTY

C01-M / C02-M LED CURING LIGHTS MINI (with USB Charger)




- 360° Rotatable tips
- Automatic compensation circuit for constant power output
- Fully covered membrane panel for easy cleaning and infection control
- Light weight: approx. 85 gms
- Dual modes: FULL / SOFT
- 12 Operation times
- Fully charged battery provides over 100 x 10 second curing cycles

C01-M / C02-M SPECIFICATIONS		
Model	C01-M (F.O. Light Guide)	C02-M (Direct Light Source)
LIGHT TIP	8 mm (dia.)	11 mm (dia.)
LIGHT WAVELENGTH	440~480 nm	440~480 nm
MODE	FULL	1,200 mW/cm ² for 5, 10, 15, 20, 25, 30 seconds±10%
	SOFT	1 st five seconds from 0-900 mW/cm ² next 5, 10, 15, 20, 25 seconds 900 mW/cm ² ±10%
WEIGHT	86 gms	84 gms
LENGTH	230 mm	220 mm

C01-M / C02-M CURING LIGHT DEPTH OF CURE							
MODE	MODEL	C01-M (F.O. Light Guide) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-M (Direct Light Source) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
FULL		10 seconds	4.3 mm	2.15 mm	10 seconds	5.1 mm	2.55 mm
SOFT		15 seconds	4.6 mm	2.30 mm	15 seconds	4.8 mm	2.40 mm

REMARKS: Above data is based on test results using an A2 shade composite. The light tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two (50% of the total composite depth cured). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is recommended that you should test your resins and composites performance to determine the curing times required.



C10 LED & HALOGEN CURING LIGHT METER



- Simple, compact design for easy use
- Measures intensity up to 2500 mw/cm²
- Size: 110 mm x 68 mm x 33 mm


C10 Light Meter



C01-M Curing Light (F.O. Light Guide)
 C02-M Curing Light (Direct Light Head)
 C01-M-1 F.O. Light Guide 8 mm
 C02-M-1 Direct Light Head 11 mm
 C01/2-M-1 Replacement Battery
 228 Light Sleeves *Mini* 500 pcs/box


24 MONTH WARRANTY

C01 LED CURING LIGHTS with fibre optic light guides




C01 SPECIFICATIONS			
Model	C01-C (classic)	C01-S (super power)	C01-D (dual range)
F.O. LIGHT GUIDE	8 mm (dia.)	11 mm (dia.)	
LIGHT WAVELENGTH	440~480 nm	440~480 nm	390~430 nm 440~480 nm
APPLICATION	General	Ortho / General	General / Ortho
MODE	TACK	N/A	900 mW/cm ² for 1 second±10%
	TURBO	1,500 mW/cm ² for 5 seconds ±10%	2,000 mW/cm ² for 4 seconds±10%
	3' + 3'	N/A	2,000 mW/cm ² for 3 + 3 seconds±10% (with 1 second blink in between)
	6' + 6'	N/A	1,500 mW/cm ² for 6 + 6 seconds±10% (with 1 second blink in between)
	FULL	1,200 mW/cm ² for 5, 10, 15, 20 seconds±10%	
	RAMP	1 st five seconds from 0-1,200 mW/cm ² next 5, 10, 15, 20 seconds 1,200 mW/cm ² ±10%	
	SOFT CURE	1 st five seconds from 0-900 mW/cm ² next 5, 10, 15, 20, 25 seconds 900 mW/cm ² ±10%	
WEIGHT	205 gms	200 gms	200 gms
LENGTH	260 mm	265 mm	265 mm

CURING LIGHT MODE APPLICATIONS	
TACK	Tack cure, a very short initial curing for temporary processing before final cure. For example, veneer bonding or temporary cement for in-position and easy removal of excess material.
TURBO	Super fast curing for resin and composite under 2 mm thickness. Also recommended for ceramic brackets, veneer bonding etc.
3' + 3'	3 Seconds Mesial + 3 Seconds Distal for metal brackets resin bonding. (with 1 second blink in between)
6' + 6'	6 Seconds Mesial + 6 Seconds Distal for metal brackets with bonding resin requiring more curing time, and buccal tubes bonding. (with 1 second blink in between)
FULL	For general applications.
RAMP	With first 0-5 seconds ramping for general application.
SOFT CURE	With first 0-5 seconds ramping for more gentle processing, designed for the patients sensitive to the heat.



C01-C Curing Light (Classic) with 8 mm Light Guide
 C01-S Curing Light (Super Power) with 11 mm Light Guide
 C01-D Curing Light (Dual Range) with 11 mm Light Guide
 C01-C-1 8 mm Fibre Optic Light Guide
 C01-S/D-1 11 mm Fibre Optic Light Guide
 C01-S/D-4 Tack Heat Protective Rings, 6 pcs/pack
 C01/2-1 Replacement Battery
 123S Light Sleeves Small (for C01-C) 500 pcs/box
 123L Light Sleeves Large (for C01-S/C01-D) 500 pcs/box

24 MONTH WARRANTY



C01-C / C01-S / C01-D CURING LIGHT DEPTH OF CURE									
MODEL	C01-C (classic) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C01-S (super power) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C01-D (dual range) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
TURBO	5 seconds	4.2 mm	2.10 mm	4 seconds	4.6 mm	2.30 mm	4 seconds	4.3 mm	2.15 mm
FULL	10 seconds	4.3 mm	2.15 mm	10 seconds	5.1 mm	2.55 mm	10 seconds	4.4 mm	2.20 mm
RAMP	10 seconds	4.1 mm	2.05 mm	10 seconds	5.0 mm	2.50 mm	10 seconds	4.2 mm	2.10 mm
SOFT CURE	15 seconds	4.6 mm	2.30 mm	15 seconds	4.8 mm	2.40 mm	15 seconds	4.6 mm	2.30 mm

REMARKS: Above data is based on test results using an A2 shade composite. The light guide tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two (50% of the total composite depth cured). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is strongly recommended that you should test your composite and resin performance to determine the curing times required.

Data is only for reference, always follow the instructions from your resin and composite manufacturers.

- There are many variables (distance and angle between the light tip and composite, area, shade, thickness and type etc.) affecting composite curing times. It is recommended to place and cure a maximum of 2 mm for each composite placement.
- Using a barrier sleeve on the curing light will reduce 5-10% of light intensity.
- The light guide tip should be placed approximately 1-2 mm above the composite or resin and pointed perpendicular to the tooth surface.
- High intensity curing lights produce more heat. This is a direct indication of the power of the light being emitted. When using the device for more than 10 seconds, do not keep the light tip in one position. It is always recommended to use a dental syringe to blow a cold air flow on the area during long periods of curing. This will maintain a lower temperature for the composite and surrounding area, particularly when using rubber dam so as to avoid any rubber dam breakage due to heat built-up.
- Do not expose soft tissue to the light being emitted for more than 2 seconds or burns may occur.
- Avoid the light guide tip coming into contact with any soft tissue in the mouth.

10 colours available:




C02 LED CURING LIGHTS 90° right angle direct light source head




C02 SPECIFICATIONS			
Model	C02-C (classic)	C02-S (super power)	C02-D (dual range)
DIRECT LIGHT HEAD	10 mm (dia.)		
LIGHT WAVELENGTH	440~480 nm	440~480 nm	390~430 nm 440~480 nm
APPLICATION	General	Ortho / General	General / Ortho
MODE	TACK	N/A	900 mW/cm ² for 1 second±10%
	TURBO	1,500 mW/cm ² for 5 seconds ±10%	2,500 mW/cm ² for 3 seconds±10%
	3' + 3'	N/A	2,000 mW/cm ² for 3 + 3 seconds±10% (with 1 second blink in between)
	6' + 6'	N/A	1,500 mW/cm ² for 6 + 6 seconds±10% (with 1 second blink in between)
	FULL	1,200 mW/cm ² for 5, 10, 15, 20 seconds±10%	
	RAMP	1 st five seconds from 0-1,200 mW/cm ² next 5, 10, 15, 20 seconds 1,200 mW/cm ² ±10%	
SOFT CURE	1 st five seconds from 0-900 mW/cm ² next 5, 10, 15, 20, 25 seconds 900 mW/cm ² ±10%		
WEIGHT	185 gms	185 gms	185 gms
LENGTH	255 mm	255 mm	255 mm

CURING LIGHT MODE APPLICATIONS	
TACK	Tack cure, a very short initial curing for temporary processing before final cure. For example, veneer bonding or temporary cement for in-position and easy removal of excess material.
TURBO	Super fast curing for resin and composite under 2 mm thickness. Also recommended for ceramic brackets, veneer bonding etc.
3' + 3'	3 Seconds Mesial + 3 Seconds Distal for metal brackets resin bonding. (with 1 second blink in between)
6' + 6'	6 Seconds Mesial + 6 Seconds Distal for metal brackets with bonding resin requiring more curing time, and buccal tubes bonding. (with 1 second blink in between)
FULL	For general applications.
RAMP	With first 0-5 seconds ramping for general application.
SOFT CURE	With first 0-5 seconds ramping for more gentle processing, designed for the patients sensitive to the heat.



C02-C (Classic) Curing Light with 10 mm Light Head
C02-S (Super Power) Curing Light with 10 mm Light Head
C02-D (Dual Range) Curing Light with 10 mm Light Head
C02-C-1 (Classic) Direct Light Source Head
C02-S-1 (Super Power) Direct Light Source Head
C02-D-1 (Dual Range) Direct Light Source Head
C01/2-1 Replacement Battery
C02-3 Heat Protect Rings, 6 pcs/pack
C02-4 Tack Heat Protective Rings, 6 pcs/pack
125 Light Sleeves (for C02-C/S/D) 500 pcs/box

24 MONTH WARRANTY



C02-C / C02-S / C02-D CURING LIGHT DEPTH OF CURE

MODEL	C02-C (classic) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-S (super power) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)	C02-D (dual range) Curing Time	Depth Measurements	ISO 4049 Depth Measurements (-50%)
TURBO	5 seconds	4.4 mm	2.20 mm	3 seconds	4.6 mm	2.30 mm	3 seconds	4.5 mm	2.25 mm
FULL	10 seconds	5.1 mm	2.55 mm	10 seconds	5.3 mm	2.65 mm	10 seconds	5.2 mm	2.60 mm
RAMP	10 seconds	4.8 mm	2.40 mm	10 seconds	5.3 mm	2.65 mm	10 seconds	4.9 mm	2.45 mm
SOFT CURE	15 seconds	5.0 mm	2.50 mm	15 seconds	5.4 mm	2.70 mm	15 seconds	5.1 mm	2.55 mm

REMARKS: Above data is based on test results using an A2 shade composite. The light tip was placed above and pointed perpendicular to the composite. Tests were done as per ISO 4049 Standard. The samples were measured and the data was divided by two (50% of the total composite depth cured). ISO 4049 Standard requires that the results for class II restorative materials should be at least 1.5 mm for non-opaque shades and 1.0 mm for opaque shades. All resins and composites have different curing times. It is strongly recommended that you should test your composite and resin performance to determine the curing times required.


Data is only for reference, always follow the instructions from your resin and composite manufacturers.

- There are many variables (distance and angle between the light tip and composite, area, shade, thickness and type etc.) affecting composite curing times. It is recommended to place and cure a maximum of 2 mm for each composite placement.
- Using a barrier sleeve on the curing light will reduce 5-10% of light intensity.
- The light tip should be placed approximately 1-2 mm above the composite or resin and pointed perpendicular to the tooth surface.
- High intensity curing lights produce more heat. This is a direct indication of the power of the light being emitted. When using the device for more than 10 seconds, do not keep the light tip in one position. It is always recommended to use a dental syringe to blow a cold air flow on the area during long periods of curing. This will maintain a lower temperature for the composite and surrounding area, particularly when using rubber dam so as to avoid any rubber dam breakage due to heat build-up.
- Do not expose soft tissue to the light being emitted for more than 2 seconds or burns may occur.
- Avoid the light tip coming into contact with any soft tissue in the mouth. As an option, you may cover the entire curing light with a disposable sleeve and then place the Heat Protective Ring (included inside the packing box) onto the light tip covered with sleeve as shown below. Ensure the Heat Protective Ring is securely positioned on the curing light before starting the cure.

10 colours available:



Prophy Angles



DISPOSABLE 60° TURN RECIPROCATING LATEX FREE PROPHY ANGLES

NO SPLATTER - The centrifugal force generated by a normal 360° turn prophy angle results in splattering and mess. Our **NEW 60° turn reciprocating** prophy angles eliminates potential splattering of prophy paste, blood and saliva, keeping the prophy paste in the cup where it belongs.

NO HEAT - Our **NEW 60° turn reciprocating** prophy angles build up much less heat during use, assuring a smooth, cool treatment for your patient.

EFFECTIVE - Due to less heat generated during use, if necessary more pressure may be applied to remove stubborn stains.

FASTER - Thanks to our innovatively designed **NEW 60° turn reciprocating** prophy angles, the prophy cup or brush can be used continuously on the tooth surface, saving valuable chair time.

SAFER - The rotational motion of traditional prophy cups and brushes can cause laceration and bleeding of the gums. Our **NEW 60° turn reciprocating** prophy angles provide a much safer and pain-free prophylaxis treatment.

WARNING: USE OF THESE 60° TURN RECIPROCATING PROPHY ANGLES IS EXTREMELY ADDICTIVE
Once you've tried them, you will never want to use any other prophy angles again

DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES STANDARD CUPS LATEX FREE TURBINE BLADES



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	CT3399-60WT	CT3399-60LB	CT3399-60LV	CT3399-60OR
Regular Assorted Colours Order No.	CT3399-60RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	CT3399-60PK	CT3399-60LM	CT3399-60GR	CT3399-60DB
Soft Assorted Colours Order No.	CT3399-60SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	CT3399-60RD	CT3399-60YL
X-Soft Assorted Colours Order No.	CT3399-60XA Assorted: Red, Yellow	

DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES SHORT CUPS LATEX FREE TURBINE BLADES



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	CTS3399-60WT	CTS3399-60LB	CTS3399-60LV	CTS3399-60OR
Regular Assorted Colours Order No.	CTS3399-60RA Assorted: White, L.Blue, Lavender, Orange			



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	CTS3399-60PK	CTS3399-60LM	CTS3399-60GR	CTS3399-60DB
Soft Assorted Colours Order No.	CTS3399-60SA Assorted: Pink, Lime, Green, D.Blue			



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	CTS3399-60RD	CTS3399-60YL
X-Soft Assorted Colours Order No.	CTS3399-60XA Assorted: Red, Yellow	

DISPOSABLE 60° TURN RECIPROCATING PROPHY ANGLES LATEX FREE BRUSHES REGULAR / SOFT / X-SOFT




Regular Type	Flat	Tapered	Mini Cup
Regular White Colour Order No.	C3399-60-R261F	C3399-60-R261T	C3399-60-R269
Regular Assorted Colours Order No.	C3399-60-R261FA C3399-60-R261TA C3399-60-R269A Assorted: White, L.Blue, Lavender, Orange		



Soft Type	Flat	Tapered	Mini Cup
Soft Clear Colour Order No.	C3399-60-S261F	C3399-60-S261T	C3399-60-S269
Soft Assorted Colours Order No.	C3399-60-S261FA C3399-60-S261TA C3399-60-S269A Assorted: Pink, Lime, Green, D.Blue		



X-Soft Type	Flat	Tapered	Mini Cup
X-Soft Yellow Colour Order No.	C3399-60-X261F	C3399-60-X261T	C3399-60-X269
X-Soft Assorted Colours Order No.	C3399-60-X261FA C3399-60-X261TA C3399-60-X269A Assorted: Yellow, Red		



STANDARD STRAIGHT HANDPIECE 1:1

L10-1 Standard Straight Handpiece 1:1

DISPOSABLE PROPHY ANGLES STANDARD CUPS LATEX FREE RIBBED AND WEBBED



Type	90° Right Angle		105° Contra Angle	
Cup Colours	Regular Lavender	Soft Green	Regular Lavender	Soft Green
Order No.	3399-090RL	3399-090SG	3399-105RL	3399-105SG

DISPOSABLE PROPHY ANGLES BRUSHES LATEX FREE FLAT AND TAPERED



Type	90° Right Angle		105° Contra Angle	
Brushes	Regular Flat	Regular Tapered	Regular Flat	Regular Tapered
Order No.	3399-090F	3399-090T	3399-105F	3399-105T

DISPOSABLE CLEAR PROPHY ANGLES STANDARD CUPS LATEX FREE RIBBED AND WEBBED



Type	90° Right Angle		90° Right Angle		Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink	Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	C3399-090RB	C3399-090RL	C3399-090SG	C3399-090SP	Order No.	C3399-105RB	C3399-105RL	C3399-105SG	C3399-105SP

DISPOSABLE CLEAR PROPHY ANGLES SHORT CUPS LATEX FREE RIBBED AND WEBBED



Type	90° Right Angle		90° Right Angle		Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink	Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	CS3399-090RB	CS3399-090RL	CS3399-090SG	CS3399-090SP	Order No.	CS3399-105RB	CS3399-105RL	CS3399-105SG	CS3399-105SP

DISPOSABLE CLEAR PROPHY ANGLES STANDARD CUPS LATEX FREE TURBINE BLADES



Type	90° Right Angle		90° Right Angle		Type	105° Contra Angle		105° Contra Angle	
Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink	Cup Colours	Regular Light Blue	Regular Lavender	Soft Green	Soft Pink
Order No.	CT3399-090RB	CT3399-090RL	CT3399-090SG	CT3399-090SP	Order No.	CT3399-105RB	CT3399-105RL	CT3399-105SG	CT3399-105SP

Prophy Cups and Points

METAL SHANKS CUPS (8.5 mm)



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	M222WT	M222LB	M222LV	M222OR
Regular Assorted Colours Order No.	M222RA Assorted: White, L.Blue, Lavender, Orange			

METAL SHANKS CUPS (8.5 mm)



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	M222PK	M222LM	M222GR	M222DB
Soft Assorted Colours Order No.	M222SA Assorted: Pink, Lime, Green, D.Blue			

METAL SHANKS CUPS (8.5 mm)



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	M222RD	M222YL
X-Soft Assorted Colours Order No.	M222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA METAL SHANKS RIBBED AND WEBBED

PLASTIC SHANKS CUPS (8.5 mm)



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	P222WT	P222LB	P222LV	P222OR
Regular Assorted Colours Order No.	P222RA Assorted: White, L.Blue, Lavender, Orange			

PLASTIC SHANKS CUPS (8.5 mm)



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	P222PK	P222LM	P222GR	P222DB
Soft Assorted Colours Order No.	P222SA Assorted: Pink, Lime, Green, D.Blue			

PLASTIC SHANKS CUPS (8.5 mm)



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	P222RD	P222YL
X-Soft Assorted Colours Order No.	P222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA PLASTIC SHANKS RIBBED AND WEBBED

METAL SHANKS SHORT CUPS (7.0 mm)



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	S222WT	S222LB	S222LV	S222OR
Regular Assorted Colours Order No.	S222RA Assorted: White, L.Blue, Lavender, Orange			

METAL SHANKS SHORT CUPS (7.0 mm)



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	S222PK	S222LM	S222GR	S222DB
Soft Assorted Colours Order No.	S222SA Assorted: Pink, Lime, Green, D.Blue			

METAL SHANKS SHORT CUPS (7.0 mm)



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	S222RD	S222YL
X-Soft Assorted Colours Order No.	S222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS CA METAL SHANKS RIBBED AND WEBBED

SHORT SHANKS SHORT CUPS (7.0 mm)



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	SS-222WT	SS-222LB	SS-222LV	SS-222OR
Regular Assorted Colours Order No.	SS-222RA Assorted: White, L.Blue, Lavender, Orange			

SHORT SHANKS SHORT CUPS (7.0 mm)



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	SS-222PK	SS-222LM	SS-222GR	SS-222DB
Soft Assorted Colours Order No.	SS-222SA Assorted: Pink, Lime, Green, D.Blue			

SHORT SHANKS SHORT CUPS (7.0 mm)



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	SS-222RD	SS-222YL
X-Soft Assorted Colours Order No.	SS-222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS CA METAL SHANKS SHORT SHANKS RIBBED AND WEBBED

[fit mini head CA handpieces only]

REGULAR POINTS STANDARD SHANKS



Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	225WT	225LB	225LV	225OR
Regular Assorted Colours Order No.	225RA Assorted: White, L.Blue, Lavender, Orange			

SOFT POINTS STANDARD SHANKS



Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	225PK	225LM	225GR	225DB
Soft Assorted Colours Order No.	225SA Assorted: Pink, Lime, Green, D.Blue			

X-SOFT POINTS STANDARD SHANKS



X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	225RD	225YL
X-Soft Assorted Colours Order No.	225XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE PROPHY POINTS CA METAL SHANKS

[for lingual and interproximal]

METAL SHANKS CUPS (8.5 mm)

REGULAR STANDARD CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-M222 WT	T-M222 LB	T-M222 LV	T-M222 OR
Regular Assorted Colours Order No.	T-M222RA Assorted: White, L.Blue, Lavender, Orange			

METAL SHANKS CUPS (8.5 mm)

SOFT STANDARD CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-M222 PK	T-M222 LM	T-M222 GR	T-M222 DB
Soft Assorted Colours Order No.	T-M222SA Assorted: Pink, Lime, Green, D.Blue			

METAL SHANKS CUPS (8.5 mm)

X-SOFT STANDARD CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

8.5mm
15mm

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-M222 RD	T-M222 YL
X-Soft Assorted Colours Order No.	T-M222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA METAL SHANKS TURBINE BLADES

PLASTIC SHANKS CUPS (8.5 mm)

REGULAR STANDARD CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-P222 WT	T-P222 LB	T-P222 LV	T-P222 OR
Regular Assorted Colours Order No.	T-P222RA Assorted: White, L.Blue, Lavender, Orange			

PLASTIC SHANKS CUPS (8.5 mm)

SOFT STANDARD CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-P222 PK	T-P222 LM	T-P222 GR	T-P222 DB
Soft Assorted Colours Order No.	T-P222SA Assorted: Pink, Lime, Green, D.Blue			

PLASTIC SHANKS CUPS (8.5 mm)

X-SOFT STANDARD CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

8.5mm
15mm

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-P222 RD	T-P222 YL
X-Soft Assorted Colours Order No.	T-P222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS CA PLASTIC SHANKS TURBINE BLADES

METAL SHANKS SHORT CUPS (7.0 mm)

REGULAR SHORT CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-S222 WT	T-S222 LB	T-S222 LV	T-S222 OR
Regular Assorted Colours Order No.	T-S222RA Assorted: White, L.Blue, Lavender, Orange			

METAL SHANKS SHORT CUPS (7.0 mm)

SOFT SHORT CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-S222 PK	T-S222 LM	T-S222 GR	T-S222 DB
Soft Assorted Colours Order No.	T-S222SA Assorted: Pink, Lime, Green, D.Blue			

METAL SHANKS SHORT CUPS (7.0 mm)

X-SOFT SHORT CUPS
TURBINE BLADES

1.500 - 3.000 r.p.m

7 mm
15mm

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-S222 RD	T-S222 YL
X-Soft Assorted Colours Order No.	T-S222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS CA METAL SHANKS TURBINE BLADES

SHORT SHANKS SHORT CUPS (7.0 mm)

REGULAR SHORT CUPS
SHORT SHANKS TURBINE BLADES

1.500 - 3.000 r.p.m

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	SS-T222WT	SS-T222LB	SS-T222LV	SS-T222OR
Regular Assorted Colours Order No.	SS-T222RA Assorted: White, L.Blue, Lavender, Orange			

SHORT SHANKS SHORT CUPS (7.0 mm)

SOFT SHORT CUPS
SHORT SHANKS TURBINE BLADES

1.500 - 3.000 r.p.m

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	SS-T222PK	SS-T222LM	SS-T222GR	SS-T222DB
Soft Assorted Colours Order No.	SS-T222SA Assorted: Pink, Lime, Green, D.Blue			

SHORT SHANKS SHORT CUPS (7.0 mm)

X-SOFT SHORT CUPS
SHORT SHANKS TURBINE BLADES

1.500 - 3.000 r.p.m

7 mm
SHORT SHANKS 12 mm

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	SS-T222RD	SS-T222YL
X-Soft Assorted Colours Order No.	SS-T222XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS CA METAL SHANKS SHORT SHANKS TURBINE BLADES [fit mini head CA handpieces only]

REGULAR POINTS SCREW-IN

1.500 - 3.000 r.p.m

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	226WT	226LB	226LV	226OR
Regular Assorted Colours Order No.	226RA Assorted: White, L.Blue, Lavender, Orange			

SOFT POINTS SCREW-IN

1.500 - 3.000 r.p.m

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	226PK	226LM	226GR	226DB
Soft Assorted Colours Order No.	226SA Assorted: Pink, Lime, Green, D.Blue			

X-SOFT POINTS SCREW-IN

1.500 - 3.000 r.p.m

9.5 mm

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	226RD	226YL
X-Soft Assorted Colours Order No.	226XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE PROPHY POINTS SCREW-IN [for lingual and interproximal]

Prophy Brushes and Cups

DISPOSABLE LATEX FREE PROPHY BRUSHES CA REGULAR / SOFT / X-SOFT



Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Regular White Colour Order No.	R-262F	R-262T	R-265	R-M268N	R-269N
Regular Assorted Colours Order No.	R-262FA R-262TA R-265A R-268A R-269A				
	Assorted: White, L.Blue, Lavender, Orange				



Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Soft Clear Colour Order No.	S-262F	S-262T	S-265	S-268	S-269
Soft Assorted Colours Order No.	S-262FA S-262TA S-265A S-268A S-269A				
	Assorted: Pink, Lime, Green, D.Blue				



X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
X-Soft Yellow Colour Order No.	X-262F	X-262T	X-265	X-268	X-269
X-Soft Assorted Colours Order No.	X-262FA X-262TA X-265A X-268A X-269A				
	Assorted: Red, Yellow				

DISPOSABLE LATEX FREE SHORT SHANKS PROPHY BRUSHES CA [fit mini head CA handpieces only]



Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Regular White Colour Order No.	SS-R262F	SS-R262T	SS-R265	SS-R268	SS-R269
Regular Assorted Colours Order No.	SS-R262FA SS-R262TA SS-R265A SS-R268A SS-R269A				
	Assorted: White, L.Blue, Lavender, Orange				



Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Soft Clear Colour Order No.	SS-S262F	SS-S262T	SS-S265	SS-S268	SS-S269
Soft Assorted Colours Order No.	SS-S262FA SS-S262TA SS-S265A SS-S268A SS-S269A				
	Assorted: Pink, Lime, Green, D.Blue				



X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
X-Soft Yellow Colour Order No.	SS-X262F	SS-X262T	SS-X265	SS-X268	SS-X269
X-Soft Assorted Colours Order No.	SS-X262FA SS-X262TA SS-X265A SS-X268A SS-X269A				
	Assorted: Red, Yellow				

SNAP-ON CUPS (8.0 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

REGULAR CUPS RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	220WT	220LB	220LV	220OR
Regular Assorted Colours Order No.	220RA Assorted: White, L.Blue, Lavender, Orange			

SNAP-ON CUPS (8.0 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

SOFT CUPS RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	220PK	220LM	220GR	220DB
Soft Assorted Colours Order No.	220SA Assorted: Pink, Lime, Green, D.Blue			

SNAP-ON CUPS (8.0 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

X-SOFT CUPS RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	220RD	220YL
X-Soft Assorted Colours Order No.	220XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE PROPHY CUPS SNAP-ON RIBBED AND WEBBED

SCREW-IN CUPS (8.5 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

REGULAR STANDARD CUPS RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	221WT	221LB	221LV	221OR
Regular Assorted Colours Order No.	221RA Assorted: White, L.Blue, Lavender, Orange			

SCREW-IN CUPS (8.5 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

SOFT STANDARD CUPS RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	221PK	221LM	221GR	221DB
Soft Assorted Colours Order No.	221SA Assorted: Pink, Lime, Green, D.Blue			

SCREW-IN CUPS (8.5 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

X-SOFT STANDARD CUPS RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	221RD	221YL
X-Soft Assorted Colours Order No.	221XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE STANDARD PROPHY CUPS SCREW-IN RIBBED AND WEBBED

SCREW-IN SHORT CUPS (7.0 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

REGULAR SHORT CUPS RIBBED AND WEBBED

Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	S221WT	S221LB	S221LV	S221OR
Regular Assorted Colours Order No.	S221RA Assorted: White, L.Blue, Lavender, Orange			

SCREW-IN SHORT CUPS (7.0 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

SOFT SHORT CUPS RIBBED AND WEBBED

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	S221PK	S221LM	S221GR	S221DB
Soft Assorted Colours Order No.	S221SA Assorted: Pink, Lime, Green, D.Blue			

SCREW-IN SHORT CUPS (7.0 mm)

RIBBED AND WEBBED

1,500 - 3,000 r.p.m.

X-SOFT SHORT CUPS RIBBED AND WEBBED

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	S221RD	S221YL
X-Soft Assorted Colours Order No.	S221XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS SCREW-IN RIBBED AND WEBBED

DISPOSABLE LATEX FREE PROPHY BRUSHES SNAP-ON REGULAR / SOFT / X-SOFT



REGULAR WHITE
(nylon brush dia. 0.20 mm)

1,500 - 3,000 r.p.m.



SOFT CLEAR
(nylon brush dia. 0.15 mm)

1,500 - 3,000 r.p.m.



X-SOFT YELLOW
(nylon brush dia. 0.10 mm)

1,500 - 3,000 r.p.m.

Regular Type	Flat	Tapered	Soft Type	Flat	Tapered	X-Soft Type	Flat	Tapered
Regular White Brushes Order No.	R-260F	R-260T	Soft Clear Brushes Order No.	S-260F	S-260T	X-Soft Yellow Brushes Order No.	X-260F	X-260T

DISPOSABLE LATEX FREE PROPHY BRUSHES SCREW-IN REGULAR / SOFT / X-SOFT



REGULAR WHITE (nylon brush dia. 0.20 mm)

1,500 - 3,000 r.p.m.



SOFT CLEAR (nylon brush dia. 0.15 mm)


1,500 - 3,000 r.p.m.



X-SOFT YELLOW (nylon brush dia. 0.10 mm)

1,500 - 3,000 r.p.m.


Regular Type	Flat	Tapered	Micro	Junior Cup	Mini Cup	Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup	X-Soft Type	Flat	Tapered	Micro	Junior Cup	Mini Cup
Regular White Colour Order No.	R-261F	R-261T	R-265-1	R-268-1	R-269-1	Soft Clear Colour Order No.	S-261F	S-261T	S-265-1	S-268-1	S-269-1	X-Soft Yellow Colour Order No.	X-261F	X-261T	X-265-1	X-268-1	X-269-1
Regular Assorted Colours Order No.	R-261FA R-261TA R-265-1A R-268-1A R-269-1A					S-261FA S-261TA S-265-1A S-268-1A S-269-1A					X-261FA X-261TA X-265-1A X-268-1A X-269-1A						
	Assorted: White, L.Blue, Lavender, Orange					Assorted: Pink, Lime, Green, D.Blue					Assorted: Red, Yellow						



REGULAR CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.


Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-220WT	T-220LB	T-220LV	T-220OR
Regular Assorted Colours Order No.	T-220RA Assorted: White, L.Blue, Lavender, Orange			



SOFT CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-220PK	T-220LM	T-220GR	T-220DB
Soft Assorted Colours Order No.	T-220SA Assorted: Pink, Lime, Green, D.Blue			




X-SOFT CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-220RD	T-220YL
X-Soft Assorted Colours Order No.	T-220XA Assorted: Red, Yellow	


DISPOSABLE LATEX FREE PROPHY CUPS SNAP-ON TURBINE BLADES



REGULAR STANDARD CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.


Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-221WT	T-221LB	T-221LV	T-221OR
Regular Assorted Colours Order No.	T-221RA Assorted: White, L.Blue, Lavender, Orange			



SOFT STANDARD CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-221PK	T-221LM	T-221GR	T-221DB
Soft Assorted Colours Order No.	T-221SA Assorted: Pink, Lime, Green, D.Blue			




X-SOFT STANDARD CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-221RD	T-221YL
X-Soft Assorted Colours Order No.	T-221XA Assorted: Red, Yellow	


DISPOSABLE LATEX FREE STANDARD PROPHY CUPS SCREW-IN TURBINE BLADES



REGULAR SHORT CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.


Regular Cup Colours	White	Light Blue	Lavender	Orange
Regular Single Colour Order No.	T-S221WT	T-S221LB	T-S221LV	T-S221OR
Regular Assorted Colours Order No.	T-S221RA Assorted: White, L.Blue, Lavender, Orange			



SOFT SHORT CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

Soft Cup Colours	Pink	Lime	Green	Dark Blue
Soft Single Colour Order No.	T-S221PK	T-S221LM	T-S221GR	T-S221DB
Soft Assorted Colours Order No.	T-S221SA Assorted: Pink, Lime, Green, D.Blue			




X-SOFT SHORT CUPS TURBINE BLADES

1,500 - 3,000 r.p.m.

X-Soft Cup Colours	Red	Yellow
X-Soft Single Colour Order No.	T-S221RD	T-S221YL
X-Soft Assorted Colours Order No.	T-S221XA Assorted: Red, Yellow	

DISPOSABLE LATEX FREE SHORT PROPHY CUPS SCREW-IN TURBINE BLADES



Greenline air-water syringe tips

Reusable products are a great idea. Unfortunately, trying to clean the inside of a reusable metal air-water syringe tip is just not possible and if you can't clean it, you can't guarantee the effectiveness of its sterilisation. This means that using a disposable plastic air-water syringe tip is the only way to guarantee avoiding cross contamination.

- Directly interchangeable with American type syringe tips (long), or European syringes which have been converted to use American type syringe tips (short)
- Made from 100% recycled plastic
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air
- Cardboard packaging which is easily recyclable

Greenline Air-water Syringe Tips
G-018 (American Type) 200 pcs/box
G-028 (European Type) 200 pcs/box

1 Our Greenline air-water syringe tips are made from 100% recycled plastic, using waste materials left over from our other manufacturing processes.

2 The Greenline air-water syringe tips are packaged in 100% recyclable cardboard boxes, instead of large plastic bags, helping to reduce the amount of plastic packaging in your surgery.



Greenline bamboo cups

- Leakproof construction, non-slip design
- Dispenses easily from cup dispenser due to their rolled lip
- Conveniently designed 170ml (6oz) size

Looking to cut down on plastic use? Switching to bamboo cups is a simple and easy change that can have a big impact on reducing plastic waste in your surgery.

Here are some reasons to make the switch to using bamboo:

Bamboo Cups 170 ml (6 oz), 1000 pcs/box
6800BB Blank
6800BH Healthy Teeth (4 assorted designs)

1 Bamboo is the world's quickest growing plant and provides a sustainable source of material. When bamboo is cut down to make paper, it is able to quickly recover to be harvested again.

2 Bamboo absorbs a large amount of CO² as it grows, which helps to offset human carbon emissions.

3 Our Greenline outer carton is 100% recyclable.



HOLLOW HANDLE micro applicators

In an effort to limit the environmental effect of one of the most used items in a dental surgery we came up with some great ideas:



- Bends into any angle for the placement of small controlled amounts of material in limited access areas.
- Economical - saving up to 80% by the reduction of excess materials used.

Hollow Handle Micro Applicators
 (Length 100 mm, dia 3 mm handles)
 HH-901 Regular (dia. 2.5 mm) Light Green
 HH-902 Fine (dia. 2 mm) Medium Green
 HH-903 Ultra Fine (dia. 1.5 mm) Dark Green

1 45% less plastic used in each hollow handle applicator than our standard micro applicator, due to the hollow design of the handle.


2 Instead of using new plastic, we use 100% recycled plastic for the brush tip and 50% recycled plastic for the hollow handle. This recycled plastic is made using waste materials left over from our other manufacturing processes.

3 No more plastic tube that can use up to a third of the amount of plastic as the actual product - just a 100% recyclable box.



mini micro applicators

In an effort to limit the environmental effect of one of the most used items in a dental surgery we came up with some great ideas:



- Bends into any angle for the placement of small controlled amounts of material in limited access areas.
- Economical - saving up to 80% by the reduction of excess materials used.

Mini Micro Applicators
 (Length 85 mm, dia 2.5 mm handles)
 G-9901 Regular (dia. 2.5 mm) Light Green
 G-9902 Fine (dia. 2 mm) Medium Green
 G-9903 Ultra Fine (dia. 1.5 mm) Dark Green

1 We designed the handle to be a smaller, more practical size that reduces the amount of plastic used by 50% in comparison to our standard micro applicators.

2 Instead of new plastic, we used 100% recycled plastic. This recycled plastic is made using waste materials from our other manufacturing processes.

3 No more plastic tube that can use up to a third of the amount of plastic as the actual product - just a 100% recyclable box.



micro applicator tips

Our most environmentally friendly micro applicators:

- For the placement of small controlled amounts of material in limited access areas.
- Economical - saving up to 80% by the reduction of excess materials used.

Greenline Micro Applicator Tips
 G-921 Regular (dia. 2.5 mm) Light Green
 G-922 Fine (dia. 2 mm) Medium Green
Autoclavable Micro Applicator Tip Handle
 7502W White
 7502T Transparent

1 By using a micro applicator tip with a reusable handle, you can reduce the amount the plastic that you would normally use by 80% when compared to a standard micro applicator.

2 Using an autoclavable handle saves plastic. Clean it, autoclave it and re-use it. Simple and effective.

3 Our micro applicator tips are packaged in a 100% recyclable cardboard box. Recycling can prevent pollution by reducing the need to collect new raw materials.

4 Instead of using new plastic for the micro applicator tips, we use 100% recycled plastic. The recycled plastic is made using waste materials from our other manufacturing processes.

DISPOSABLE AIR-WATER SYRINGE TIPS (SHORT) EURO TYPE




- Directly interchangeable with American type short syringe tips (syringes which have been converted to use American type short syringe tips)
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air

Syringe Tips Short, 250 pcs/pack
 028S Syringe Tips *Short Spectrum*
 028W Syringe Tips *Short White*

DISPOSABLE AIR-WATER 'e' SYRINGE TIPS WITHOUT CORE

ready to use, no converter required for American style syringe



- Directly interchangeable with most American style metal air/water syringe tips
- Economically designed without separate water core, in exchange for a split second of moisture when switching from water to air
- Rounded tip end provides maximum patient comfort

Disposable 'e' Air-Water Syringe Tips without core, 250 pcs/pack
 0120S Opaque Colour Tips
 012W Opaque White Tips

DISPOSABLE AIR-WATER 'a' SYRINGE TIPS ready to use, no converter required fo American style syringe



- Directly interchangeable with most American style metal air/water syringe tips
- No sharp edges for patient comfort
- Unique quad air core design to ensure dry air

Disposable 'a' Air-Water Syringe Tips, 250 pcs/pack
 018S *Spectrum* Colour Tips/White Cores
 018C *Clear* Tips/Colour Cores
 018OS *Opaque Colour* Tips/White Cores
 018W *Opaque White* Tips/White Cores

DISPOSABLE TRANSPARENT SYRINGE TIPS



- Directly interchangeable with *Sani-Tip® type transparent tip (syringes which have been converted to use transparent tip)
- Costs less than sterilizing metal tips
- Size: Long(American type) Short(European type)
- Colour: Clear Transparent

Clear Tips, 250 pcs/pack
 011A Standard Length (American Type Syringe)
 011E Short Length (European Type Syringe)

*Sani-Tip is a registered trademark of Dentsply Sirona

DISPOSABLE FACE MASKS EAR-LOOP WITH PATTERNS FLUID RESISTANT LATEX FREE



LEVEL 3 BARRIER
INTERNATIONAL STANDARD ASTM F2100
BFE ≥99% PFE ≥99%

PERFORMANCE TYPE IIR
EUROPEAN STANDARD EN14683

- Non-irritating
- Fluid resistant, excellent filtration (BFE 99%)
- Permits air to pass through mask to prevent glasses from fogging

2288 Face Masks Ear-Loop with Patterns, 50 pcs/box

DISPOSABLE FLUID RESISTANT DENTAL BIBS WITH PATTERNS




- Two plies of facial quality tissue plus a layer of polyethylene film offers quick absorbency with a moisture barrier
- Paper specially treated for wet strength
- Embossed horizontal pattern, water contact disperses horizontally instead of "channelling" down onto patients clothing
- Size: 45 cm x 33 cm

Dental Bibs with Patterns

2200FA Flowers Assorted, 500 pcs/box
2200LA Leaves Assorted, 500 pcs/box
2200KID Kids, 500 pcs/box
2200MON Monsters, 500 pcs/box


DISPOSABLE PROPHY BIBS EXTRA LONG FLUID RESISTANT



- Ideal for prophy
- Features contour neck design for additional coverage to shoulders and arm
- Two plies of facial quality tissue plus a layer of polyethylene
- Size: W45 cm x L60 cm
- Colour: Light Blue

2208 Prophy Bibs Extra Long
100 pcs/box + Bib Holder XS


DISPOSABLE FLUID RESISTANT TRAY COVER PAPERS



- Two plies of facial quality tissue plus a layer of polyethylene film
- Colours: Blue and white barrier

2220 "B" Plastic Tray Covers 200'S (29.5 cm X 21.5 cm)
2221 "F" Plastic Tray Covers 200'S (20 cm X 13.5 cm)
2223 "ISO" Aluminium Tray Covers 200'S (27.5 cm X 17.5 cm)


DISPOSABLE NAPKIN HOLDERS



- Fully disposable
- Easy to apply
- Width: 15 mm
- Length: 44 cm extendable up to 60 cm

2209 Napkin Holders, 100 pcs/pack

DISPOSABLE PSP BARRIER ENVELOPES



- Quick tear-off opening
- Fits most phosphor plates
- Soft vinyl with rounded corners

DURR® TYPE

Mid-Opening Barrier Envelopes
 188-0 No. 0 (I.D. 37 mm x 24 mm) 100 pcs/box
 188-1 No. 1 (I.D. 42 mm x 26 mm) 100 pcs/box
 188-2 No. 2 (I.D. 43 mm x 33 mm) 100 pcs/box
 188-3 No. 3 (I.D. 56 mm x 29 mm) 50 pcs/box
 188-4 No. 4 (I.D. 78 mm x 59 mm) 50 pcs/box
Mid-Opening Barrier Envelopes Bulk
 188-2B No. 2 (I.D. 43 mm x 33 mm) 500 pcs/box

Top-Opening Barrier Envelopes

198-0 No. 0 (I.D. 24 mm x 37 mm) 100 pcs/box
 198-1 No. 1 (I.D. 26 mm x 42 mm) 100 pcs/box
 198-2 No. 2 (I.D. 33 mm x 43 mm) 100 pcs/box
Top-Opening Barrier Envelopes Bulk
 198-2B No. 2 (I.D. 33 mm x 43 mm) 500 pcs/box

DIGORA® (SOREDEX®) TYPE

Top-Opening Barrier Envelopes
 189-0 No. 0 (I.D. 24 mm x 34 mm) 100 pcs/box
 189-1 No. 1 (I.D. 26 mm x 43 mm) 100 pcs/box
 189-2 No. 2 (I.D. 33 mm x 44 mm) 100 pcs/box

KODAK® (CARESTEAM®) TYPE

Top-Opening Barrier Envelopes
 199-0 No. 0 (I.D. 31 mm x 36 mm) 100 pcs/box
 199-1 No. 1 (I.D. 33 mm x 43 mm) 100 pcs/box
 199-2 No. 2 (I.D. 40 mm x 46 mm) 100 pcs/box

Protective Paper Covers

178-0 (2 x 23 mm x 36 mm) 100 pcs/box
 178-1 (2 x 24 mm x 42 mm) 100 pcs/box
 178-2 (2 x 30 mm x 42 mm) 100 pcs/box

DISPOSABLE PC MOUSE SLEEVES



- Universal size, fit both cordless and corded optical mouse
- Soft, pliable polypropylene sleeves
- Easy to apply

PC Mouse Sleeves, 500 pcs/box
 133S 15 cm x 8.5 cm
 133L 16 cm x 10 cm

DISPOSABLE X-RAY SLEEVES



- Full x-ray unit coverage
- Disposable yet economical
- Soft, pliable polypropylene sleeve holds position until removed
- Easy application and removal

X-Ray Head Sleeves, 250 pcs/box
 108 Small 20 cm x 17 cm x 60 cm
 109 Large 35 cm x 24 cm x 75 cm



DISPOSABLE HVE TUBE SLEEVES

- Pre-cut opening for easier application
- Made from very supple poly plastic for greater tactile feel

HVE Suction Tube Sleeves, 500 pcs/box
 135S Small (11 mm) 5.5 cm x 35 cm
 135L Large (16 mm) 6 cm x 38 cm



A. STANDARD

140XS
 140S
 140M

A. Standard Sensor Sleeve, 500 pcs/box
 140XS X-Small 3.8 cm x 20 cm
 140S Small 4.2 cm x 20 cm
 140M Medium 4.5 cm x 20 cm
 140L Large 6 cm x 25 cm



B. TWO STEP

183-1
 183-2

B. Two Step Sensor Sleeves, 500 pcs/box
 183-1 Size 1 3.4 cm / 5 cm x 22 cm
 183-2 Size 2 3.8 cm / 5.5 cm x 22 cm

DISPOSABLE X-RAY SENSOR SLEEVES

- Made from very supple poly plastic for greater tactile feel
- Easy to apply



DISPOSABLE INSTRUMENT TRAY SLEEVES

- Instruments are placed directly on top of the disposable plastic sleeve
- Reduces the possibility of cross-contamination, saves cleaning time

Tray Sleeves, 500 pcs/box
 101 Standard 'B' 27 cm x 36 cm
 102 Mini 'F' 19 cm x 28 cm



DISPOSABLE ISO INSTRUMENT TRAY SLEEVES

- Instruments are placed directly on top of the disposable plastic sleeve
- Reduces the possibility of cross-contamination, saves cleaning time

ISO Metal Tray Sleeves, 500 pcs/box
 190 Standard 31 cm x 21 cm



DISPOSABLE IMPRESSION GUN SLEEVES

- Plastic sleeve covers cartridge and handle
- Provides protection in case of cartridge failure or leakage

Impression Gun Sleeves, 500 pcs/pack
 107 Universal 25 cm x 21 cm




DISPOSABLE MOUTH MIRROR SLEEVES

- Made from very supple poly plastic for greater tactile feel

Mouth Mirror Sleeves, 500 pcs/box
 140XS X-Small 3.8 cm x 20 cm

DISPOSABLE BARRIER FILM



- With 5 mm non-adhesive edge for simple placement and removal
- Residue free
- Colours: Blue and Clear

8008-B 8008-C

Barrier Film, 1200 pcs/roll
 8008-C Clear 10 cm x 15 cm
 8008-B Blue 10 cm x 15 cm

DISPOSABLE T-STYLE LIGHT HANDLE SLEEVES



- Universal sizing to fit most T-style handles
- Full handle coverage
- Disposable yet economical
- Soft, pliable polypropylene sleeve holds position until removed
- Easy application and removal

106 T-SHAPE

T-Shape Light Handle Sleeves, 500 pcs/box
 106 Universal 15 cm x 8.5/10 cm




DISPOSABLE BITE BLOCK SLEEVES

- Polypropylene sleeves
- Easy application and removal
- Fully disposable

145

Bite Block Sleeves, 5 x 200 pcs/box
 145 Universal 2.6 cm x 5.5 cm



Pen Type Curing Light Sleeves, 500 pcs/box
228 Mini Short 3.8 cm x 25 cm
123SS Small Short 5 cm x 28 cm

123S Small 5 cm x 32 cm
123L Large 6.5 cm x 35 cm
123XL X-Large 7.2 cm x 38 cm

DISPOSABLE PEN TYPE CURING LIGHT SLEEVES

- Made from very supple poly plastic for greater tactile feel
- Easy to apply

DISPOSABLE INTRAORAL CAMERA /RA CURING LIGHT SLEEVES



- Made from very supple poly plastic for greater tactile feel
- Easy to apply

Intraoral Camera Sleeves / RA Curing Light Sleeves
228 Mini Short 3.8 cm x 25 cm, 500 pcs/box
125SS Short 4.9 cm x 28 cm, 500 pcs/box
125 Standard 4.9 cm x 32 cm, 500 pcs/box



Light Guide Sleeves, 200 pcs/box
128 Small (fits 7-8 mm Guide)
129 Large (fits 9-11 mm Guide)

DISPOSABLE PEEL-N-STICK LIGHT GUIDE SLEEVES

- Soft, pliable polypropylene sleeve
- Easy to apply
- Economical



500 pcs/box
103 Pen Shape Instrument 2.4 cm x 22 cm
104S Low Speed/Universal Slim 3.8 cm x 22 cm
104L Low Speed/Universal Large 4 cm x 22 cm

DISPOSABLE PEN AND HANDPIECE SLEEVES

- Disposable sleeves for pen, handpieces, and ultrasonic scalers
- Made from soft pliable poly plastic for greater tactile feel



Syringe Sleeves, 500 pcs/box
100C Clear 26 cm x 6.5 cm
100B Blue 26 cm x 6.5 cm

DISPOSABLE SYRINGE/HVE SLEEVES

- Pre-cut opening for easier application
- Ideal for 3-way syringes, saliva ejectors and HVE valves etc

DISPOSABLE TUBING SLEEVES



118 Tubing Sleeve for straight tubing 2" (5 cm)
119 Tubing Sleeve for coiled tubing 4" (10 cm)

- For straight and coiled tubing
- Affordable protection from infectious aerosols and minimizes time spent on cleaning
- 1200 ft/roll (366 m)

DISPOSABLE FLOWABLE COMPOSITE SYRINGE SLEEVES



- Made from very supple poly plastic for greater tactile feel
- Easy to apply
- Economical

Flowable Composite Syringe Sleeves, 200 pcs/box
138 Universal 12.5 cm x 5 cm



Keyboard Sleeves
158S Standard 16 cm x 55 cm, 500 pcs/box
158M Medium 21.5 cm x 55 cm, 250 pcs/box
108 Large (20 cm + 17 cm) x 60 cm, 250 pcs/box

DISPOSABLE KEYBOARD SLEEVES

- Soft, pliable polypropylene sleeves



111
112

DISPOSABLE CHAIR SLEEVES

- Eliminates the need for separate headrest covers
- Avoids the need for harmful disinfectants and aerosols prolonging the life of dental chairs/stools


Chair Sleeves, 125 pcs/box
111 Full Chair 180 cm x 72 cm
112 Half Chair 100 cm x 65 cm

DISPOSABLE HEADREST SLEEVES



- Polypropylene sleeves
- Reduces possible cross-contamination from infectious aerosol spray

Headrest Sleeves, 500 pcs/box
113 Small 27.5 cm x 26 cm
114 Large 35 cm x 28 cm



DISPOSABLE NON-WOVEN HEADREST COVERS

- Breathable non-woven fabric
- Soft, provides maximum patient comfort
- Colours: Light Blue and White

Non Woven Headrest Covers, 500 pcs/box
780S Small 26 cm x 28 cm
780M Medium 35 cm x 28 cm

Applicators and Tips




DISPOSABLE MICRO APPLICATORS

- Economical and disposable, saving up to 80% by the reduction of excess materials used
- Bends into any angle for the placement of small controlled amounts of material in limited access areas



Available in :-
Standard Handles
100 mm x dia. 3 mm
Mini Handles
85 mm x dia. 2.5 mm



Available in :-
Standard Handles
100 mm x dia. 3 mm
Mini Handles
85 mm x dia. 2.5 mm



Available in :-
Standard Handles
100 mm x dia. 3 mm


Standard Handles
Tube Pack 400 pieces
901T Regular
902T Fine
903T Ultrafine
904T Cylinder

Mini Handles
Tube Pack 400 pieces
9901T Regular
9902T Fine
9903T Ultrafine

Standard Handles
Bulk Pack 1000 pieces
901B Regular
902B Fine
903B Ultrafine
904B Cylinder

Mini Handles
Bulk Pack 1000 pieces
9901B Regular
9902B Fine
9903B Ultrafine

Standard Handles
Dispenser Pack 400 pieces
901D Regular
902D Fine
903D Ultrafine
904D Cylinder



DISPOSABLE MICRO APPLICATOR DISPENSERS

- Small, sturdy, durable construction

Applicator Dispenser
18-900 Standard Handle
18-9900 Mini Handle

DISPOSABLE MICRO APPLICATOR TIPS




- Economical and disposable saving up to 80% by the reduction of excess materials used
- For the placement of small controlled amounts of material in limited access areas
- Can be used for etchants, bonding agents, cavity lining, fluoride varnishing, hemostatic and sealants etc

Micro Applicator Tips, 100 pcs/box
(Blue, Green, Lavender or Yellow)
921 Regular
922 Fine

Micro Applicator Tips, 500 pcs/box
(Blue, Green, Lavender or Yellow)
B921 Regular
B922 Fine


Brush Handle Autoclavable
7502W - White
7502T - Transparent

DISPOSABLE BRUSH TIPS AND HANDLES



- Perfect for placement of etchants, bonding agents, sealers and cements etc
- Available Colours: Blue, Green, Lavender or Yellow

7501 Brush Tips Regular Standard Pack, 100 pcs/box
7511 Brush Tips Mini Standard Pack, 100 pcs/box
7502W Brush Handle Autoclavable White
7502T Brush Handle Autoclavable Clear




DISPOSABLE BRUSH APPLICATORS

- Perfect for placement of etchants, bonding agents, sealers and cements etc
- Bendable to reach around corners and the tightest spaces
- Handles are colour coded to eliminate confusion between materials
- Disposable, eliminates cross-contamination
- Available colours: Blue, Green, Lavender or Yellow


Brush Applicator Single Colour, 100 pcs/pack
7500 Regular
7510 Mini

AUTOCLAVABLE HVE SUCTION TUBES

- Fits European style 16 mm or standard 11 mm suction tubing
- Single use, autoclavable up to 134°C
- Special high quality plastic, sturdy enough for cheek retraction
- Fully polished internal surface and edge for quiet operation and patient comfort
- Assortment of bright, lively colours: Lavender, Light Blue, Light Green, Lime and Yellow



Assorted Colours, 10 pcs/pack
 335L Adult Large (16 mm)
 335M Adult Medium (16 mm)
 336P Pedo (16 mm)
 005L Adult Large (16 mm)
 005M Adult Medium (16 mm)
 006P Pedo (16 mm)
 205M Adult Medium (11 mm)
 206P Pedo (11 mm)



DISPOSABLE SOFT TOUCH™ ORAL EVACUATORS (SLOTTED) VENTED LATEX FREE

- Soft Touch edge eliminates the risk of lacerating the oral tissue
- Rigid tubing strong enough for cheek retraction
- Vented to reduce unnecessary tissue grab and increase evacuation cover area
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders
- Colour: White (Blue or Lavender Tip)

PATENT PENDING

ST001 Soft Touch™ Oral Evacuators, 50 pcs/pack



DISPOSABLE VENTED (SLOTTED) ORAL EVACUATORS

- Vented to reduce unnecessary tissue grab
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders

Disposable Oral Evacuators, 100 pcs/pack
 001 White
 001A Single Colour



DISPOSABLE VENTED/NON-VENTED ORAL EVACUATORS 'S' CUT

- Vented and non-vented combination
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders


Disposable Oral Evacuators, 100 pcs/pack
 301 White
 301A Single Colour



DISPOSABLE VENTED ORAL EVACUATORS 'S' CUT

- Vented to reduce unnecessary tissue grab
- Soft rigid tubing, allows for greater patient comfort yet is strong enough for cheek retraction
- Thoroughly polished surface for smooth and quiet operation
- Fits standard 11 mm suction holders
- Colours: Blue, Yellow, Green, Lime or Lavender


Disposable Oral Evacuators, 100 pcs/pack
 308 White
 308A Single Colour



AUTOCLAVABLE SURGICAL ASPIRATOR TIPS

- Non-sterile
- Autoclavable up to 134°C
- Fits 11 mm tubing
- Colours: Standard - White
Large orifice - Blue
- Single use


Autoclavable Surgical Tips
 25 pcs/pack
 007 Mini Tip
 009 Standard Non-Vented
 008 Standard Vented
 008V Large Orifice Non-Vented
 008V Large Orifice Vented



DISPOSABLE SALIVA EJECTORS

- Tip design that eliminates suction of tissue
- Smooth fixed tip for patient safety and comfort
- Phthalate free
- Made from non-toxic PVC


Saliva Ejectors, 100 pcs/pack
 300W White Tip/Clear Tube
 300B Blue Tip/Clear Tube
 300A White Tip/Coloured Tube



AUTOCLAVABLE ASPIRATOR BRUSHES

- Autoclavable up to 134°C
- Rigid bristles for effective cleaning

020 X-Fine dia. 5 mm, 6 pcs/pack
 021 Fine dia. 8 mm, 6 pcs/pack
 022 Medium dia. 14 mm, 6 pcs/pack
 023 Large dia. 20 mm, 6 pcs/pack



AUTOCLAVABLE SUCTION TUBE ADAPTORS

- Autoclavable up to 134°C



DISPOSABLE 3-IN-1 TRAYS

- Wide occlusal plane
- Strong and tough plastic
- Specially designed walls for extra retention
- Tear resistant quality mesh
- Colour: Blue

3-in-1 Trays
 18-N401 Posterior 36 pcs/box
 18-N402 Quadrant 36 pcs/box
 18-N403 Sideless 36 pcs/box
 18-N404 Anterior 36 pcs/box
 18-N405 Full Arch 24 pcs/box

DISPOSABLE 3-IN-1 TRAYS

Stainless Steel Wire Connectors



- Stainless steel wire connector for greater access in limited space
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh
- Colour: Blue

3-in-1 Trays
 18-S401 Posterior 36 pcs/box
 18-S402 Quadrant 36 pcs/box
 18-S403 Sideless 36 pcs/box

DISPOSABLE 3-IN-1 TRAYS

Rigid Bar Connectors




- Rigid Aluminium Bar Connector
- Strong and tough
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh
- Colour: Blue


3-in-1 Trays
 18-T401 Posterior 36 pcs/box
 18-T402 Quadrant 36 pcs/box
 18-T403 Sideless 36 pcs/box

DISPOSABLE IMPRESSION TRAYS WITH RIM LOCK

for Implant Post and General Impressions



base sections can be easily removed with a bar for implant posts



rim lock reduced thickness trim lines on labial wall for easy removal in inclination post cases

- Save time and expensive laboratory charges
- Strong, rigid distortion free plastic construction
- Anatomically designed
- Smooth edges for patient comfort
- Excellent retention for impression materials whilst providing sufficient pressure for all details
- Colour: Blue

Disposable Impression Trays, 10 pcs/pack
 410 Intro (one size each)
 411-420 Single pack



AUTOCLAVABLE IMPRESSION TRAYS WITH RIM LOCK

- Autoclavable up to 134°C (275°F)
- Strong, rigid distortion free plastic construction
- Anatomically designed
- Smooth edges for patient comfort
- Excellent retention for impression materials whilst providing sufficient pressure for all details
- Colour: Light Grey

Autoclavable Impression Trays
 430 Intro 10 pcs/pack (one size each)
 431-440 Single 5 pcs/pack



DISPOSABLE 3-IN-1 TRAYS

- Wide occlusal plane
- Strong and tough plastic
- Specially designed walls for extra retention
- Tear resistant quality mesh
- Individually wrapped

3-in-1 Trays
 N401 Posterior 36 pcs/box
 N402 Quadrant 36 pcs/box
 N403 Sideless 36 pcs/box
 N404 Anterior 36 pcs/box
 N405 Full Arch 24 pcs/box

DISPOSABLE 3-IN-1 TRAYS

Wire Connectors Stainless Steel



- RIGID Stainless Steel Wire Connector
- Stainless steel wire connector for greater access in limited space
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh

3-in-1 Trays
 S401 Posterior 36 pcs/box
 S402 Quadrant 36 pcs/box
 S403 Sideless 36 pcs/box

DISPOSABLE 3-IN-1 TRAYS


Rigid Bar Connectors



- RIGID Aluminium Bar Connector
- Strong and tough
- Wide occlusal plane
- Special rigid walls designed for extra retention
- Tear resistant quality mesh

3-in-1 Trays
 T401 Posterior 36 pcs/box
 T402 Quadrant 36 pcs/box
 T403 Sideless 36 pcs/box

DISPOSABLE PLASTIC COATED METAL PLIABLE IMPRESSION TRAYS



- Disposable
- Pliable metal

AUTOCLAVABLE ALGINATE/PLASTER CURVED SPATULA

- For alginate or plaster
- Light weight, 58 gms only
- Autoclavable up to 134°C
- Stainless steel blade
- Handle colour: Blue

MS455 Curved Spatula

DISPOSABLE IMPRESSION SYRINGES

- Pre-bent for easy application
- Disposable
- 5 ml capacity

089 Impression Syringes 50 pcs/pack


FLEXIBLE MIXING BOWLS



- For effective mixing of alginates and plasters
- Sizes: 7460 Small dia. 11 cm x H8 cm
- 7470 Medium dia. 12.5 cm x H10 cm
- Available in Transparent Red, Green, Blue, Purple, Yellow and solid Blue

7460 Flexible Mixing Bowl Small
 7470 Flexible Mixing Bowl Medium

AUTOCLAVABLE ALGINATE AND PLASTER SPATULAS



- Autoclavable up to 134°C
- Fit all flexible mixing bowls
- Will not scratch the surface of bowls and clean up easily
- Flexible control in mixing with easy grip handle
- Size: Length: 190 mm

450R Regular Blue
 450XS Extra Stiff Grey

DISPOSABLE INTRAORAL SYRINGE TIPS

- Universal fitting
- Available in Yellow or White

Intraoral Tips, 100 pcs/pack
 085 Yellow for Yellow mixing tips
 086 Clear for Green and Pink mixing tips



DISPOSABLE DYNAMIC MIXING TIPS

- Fits Pentamix® type mixing machine
- Economically priced
- Colours: Red or White

7055 Dynamic Mixing Tips, 50 pcs/pack

® Pentamix is a registered trademark of 3M or 3M ESPE AG.


AUTOCLAVABLE ENDO STAND



- Autoclavable plastic stand with disposable sponge insert
- Cleans and holds endo instruments
- Silicone ring base resists sliding
- Autoclavable up to 134°C
- Colour: Blue

4833 Endo Stand Large (dia. 45 mm)
812 Endo Sponges, 50 pcs/pack

DISPOSABLE ENDO SPONGES




811
for round endo stand
size: 2" dia. x 3/8" thick
(50mm dia. x 8mm)

812
for oval and round endo stand
size: 2 1/2" dia. x 3/8" thick
(65mm dia. x 8mm)

- Disposable blue polyurethane sponges
- Colour: Blue

Endo Sponges, 50 pcs/pack
811 for Round Endo Stand
812 for Oval and Large Endo Stand


DISPOSABLE PEEL-AND-STICK ENDO CUSHIONS



802 Peel-and-Stick
Endo Cushions
50 pcs/box

- Hold, arrange and clean small instruments during endodontic or implant procedures
- Ideal for files, screws, drills and pluggers etc
- Peel-and-stick back secures to patient napkins, tray covers or anywhere handy
- Pre-arranging saves set up time
- Dense blue polyurethane foam 75 mm long

AUTOCLAVABLE ENDO CUSHION HOLDER



- Holder for used files with disposable foam cushions
- Built in endo file ruler
- Stainless steel construction
- Autoclavable up to 134°C
- Sizes: 50 mm x 50 mm x 25 mm

1180 Endo Cushion Holder
(without Cushion)
801 Square Endo Cushions
50 pcs/pack

AUTOCLAVABLE TRIANGULAR ENDO RING AND CUSHIONS



- Disposable cushions
- Colour: Blue
- Autoclavable up to 134°C

088 Autoclavable Cushion Holders
with Large and Small Rings
803 Disposable Cushions 50 pcs/box



AUTOCLAVABLE ENDO HOLDER ALUMINIUM

2800 Endo Holder


- Specially designed for ultrasonic cleaning and autoclaving of root canal instruments
- Holds 18 hand or engine endo instrument up to size 100 all lengths
- Millimeter rule for rubber stop measurements for both right or left hand use
- Adjustable cover with silicone stop usage counter for sterilization
- Size: W65 mm x D30 mm x H55 mm
- Autoclavable up to 134°C
- Colour: Blue

AUTOCLAVABLE RUBBER DAM FRAME



4866
Rubber
Dam Frame

- Proven design
- Better grip
- Enhance patient comfort
- Autoclavable up to 134°C




AUTOCLAVABLE ENDO DISPENSER PLASTIC

N055 Endo Dispenser

- Specially designed for cleaning and autoclaving of endo instruments
- Holds 18 hand or engine instruments up to size 100 (all lengths)
- Autoclavable up to 134°C
- Size: 85 mm x 52 mm x 14 mm
- Colour: Blue

AUTOCLAVABLE ENDO MEASURING BLOCK PLASTIC



4888
Endo Measuring Block


- Accurate settings
- Durable construction
- Autoclavable up to 134°C
- Size: L100 mm x D24 mm x H32 mm
- Colour: Blue



AUTOCLAVABLE ENDO RULER ALUMINIUM

- For both hand and engine instruments
- Autoclavable up to 134°C
- Gauges from 1 to 30 mm for right or left hand

4898 Endo Ruler



10005
Endo Measuring Block


AUTOCLAVABLE ENDO MEASURING BLOCK ALUMINIUM

- Durable aluminium construction
- Accurate settings
- Autoclavable up to 134°C
- Size: L150 mm x W24 mm x H33 mm
- Colour: Blue

Welcome to our new Greenline. We believe that making simple changes to the everyday products you purchase is the easiest way to help our planet. This range is our first step towards offering you the choice of more environmentally friendly products.

After scratching our heads, we've created our interdental brushes with a few 'ez' changes:


ez interdental brushes



- ISO size (colour coded)
- Latex and BPA free
- Reusable
- Cardboard box
- Non-slip, easy grip design
- Small, compact and easy to carry

Interdental Brushes, ISO Size No. 1-6
Single or Assorted Sizes, 6 pcs/box

- 1 We have designed the shape of our handle to be small and practical in size, so that it uses up to 65% less plastic than other leading brands
- 2 Instead of using new plastic, we use 100% recycled plastic for the brush handle and 50% recycled plastic for the protective cap. This recycled plastic is made using waste materials left over from other manufacturing processes.
- 3 Unlike some other brands who package interdental brushes in unnecessarily large plastic boxes and bags, we use a cardboard box which is 100% plastic free and easily recyclable



ISO Size	1	2	3	4	5	6
Reference Size	XX Fine	X Fine	Fine	Fine -Medium	Medium	Large
Made to fit hole (mm)	<0.8	0.9-1.0	1.1-1.2	1.3-1.5	1.6-1.8	>1.9
Wire Dia. (mm)	0.35	0.35	0.45	0.45	0.45	0.45
Order Code	G800-1	G800-2	G800-3	G800-4	G800-5	G800-6
Assorted Order Code	G800-A No.1-6 Assorted					



DISPOSABLE HEALTHY TEETH BIBS AND CUPS

HEALTHY TEETH DENTAL BIBS

- Fun way to promote dental health to patients
- Two plies of facial quality tissue plus a layer of polyethylene film offers quick absorbency with a positive moisture barrier
- Paper specially treated for wet strength
- Size: 33 cm x 45 cm


HEALTHY TEETH PAPER CUPS

- Leakproof construction, non-slip design
- Dispenses easily from cup dispenser due to their rolled lip
- Fun way to promote dental health to patients



Healthy Teeth
2200HT Dental Bibs, 500 pcs/pack
6800HT Paper Cups 170 ml (6 oz), 1,000 pcs/box


AUTOCLAVABLE TRAYS PLASTIC



- Autoclavable up to 134°C

650 Standard Flat 'B' Tray 34 cm x 24 cm x 2 cm
 651 Standard Divided 'B' Tray 34 cm x 24 cm x 2 cm
 652 Mini Tray Flat 'F' Tray 24 cm x 16.5 cm x 2 cm

STANDARD COLOURS
Blue
Light Blue
Green
Yellow
White
Light Beige
Lavender



AUTOCLAVABLE CABINET TRAYS PLASTIC

- Autoclavable up to 134°C
- Indispensable drawer organizer
- Scientifically designed for maximum convenience and efficiency
- Easy to clean, sterilize and disinfect
- Resists breaking, chipping and chemical attack
- Colours: Light Blue, Blue, White, Yellow, Lavender, Green and Light Beige
- Size: 9.5 cm x 19.5 cm x 2 cm

653 Cabinet Tray



658 Standard Rack 659 Mini Rack

INSTRUMENT TRAY RACKS

- Sturdy, chromium plated steel construction
- Size:
 - Standard (8 trays) W37 cm x D24 cm x H37 cm
 - Mini (6 trays) W28 cm x D17 cm x H29 cm

658 Tray Rack Standard 'B'
 659 Tray Rack Mini 'F'

AUTOCLAVABLE INSTRUMENT MAT



655 Instrument Mat

- Autoclavable
- Silicone material
- Size: L172 mm x W50 mm x H10 mm

STANDARD COLOURS	Yellow	Light Blue	Light Beige
Blue	White	Green	Lavender


AUTOCLAVABLE NEEDLE STATION STAINLESS STEEL



- Holds any dental needle sheath
- Solid stainless construction, autoclavable
- Silicone ring base resists sliding
- Size: Dia. D55 mm x H38 mm

031 Needle Station

AUTOCLAVABLE STERILIZATION CASSETTES PLASTIC



- Sturdy, durable construction
- Autoclavable

3881 Small (5 instruments)
 3882 Medium (8 instruments)

STANDARD COLOURS	Blue	White	Light Blue	Green	Lavender
------------------	------	-------	------------	-------	----------

AUTOCLAVABLE SILICONE INSTRUMENT CODE RINGS



- Autoclavable to 134°C
- Sizes: Small - for small instruments
Large - for pliers and forceps etc
- Available in single colors: Blue, Light Blue, Green, Yellow, Pink, Red, Orange, Lavender, Black and White

Silicone Instrument Rings
 657S Small, 50 pcs/box
 657L Large, 25 pcs/box

AUTOCLAVABLE BUR/ENDO ORGANIZERS



ENDO STANDS
 4810 (78 holes) for both hand/engine endo instruments 21 mm - 30 mm
 4811 (48 holes) for both hand/engine endo instruments 21 mm - 30 mm
 4812 (36 holes) for both hand/engine endo instruments 21 mm - 30 mm
 4813 12 endo vials

- Smart, flexible design for customize setups
- Autoclavable up to 134°C
- Colour: Blue

Endo Stand, W115 mm x D52 mm x H40 mm
 4810 - 78 Holes
 4811 - 48 Holes + 4 Endo Vials
 4812 - 36 Holes + Endo Clean
 4813 - 12 Endo Vials
Bur Stand, W115 mm x D52 mm x H40 mm
 4818 - 78 FG Holes
 4819 - 42 FG + 36 CA Holes

Stainless Steel Cover, W115 mm x D52 mm x H15 mm
 4830 - Sterilization Cover

Bur/Endo Box, (box only)
 4821 Organizer Box for 1 Stand W120 mm x D55 mm x H56 mm
 4823 Organizer Box for 3 Stands W160 mm x D120 mm x H56 mm
 4825 Organizer Box for 5 Stands W265 mm x D120 mm x H56 mm




JUMPING TEETH



- Cute design
- Economical giveaway for kids

8001
Jumping Teeth
24 pcs/box



090 Napkin Holder Adult (L47 cm)
090C Napkin Holder Child (L41 cm)
090XS Napkin Holder Extra Short (L37 cm)

AUTOCLAVABLE NAPKIN HOLDERS

- Autoclavable silicone up to 134°C
 - Stainless steel alligator clips
 - Non-slip design
 - Bright, attractive colours
 - Colours available
- Adult: Light Blue, Blue, Green, Red, Yellow, White, Pink and Lavender
Child: Light Blue, Lavender, Blue and Pink
XS (for neck contour bibs): Blue only

PEN / TOOTHBRUSH HOLDER



- Cute design
- 2005 Pen / Toothbrush Holder




DENTAL BIB HOLDER AND CLIP

- Cold disinfection only
- Cute design

091 Clip
092 Bib Holder, White or Red

DISPOSABLE MOUTH MIRRORS



- Disposable plastic mouth mirrors for patients
- Anti-fog front surface fiber glass mirror
- Mirror size: dia. 23 mm

800 Disposable Mouth Mirrors, 10 pcs/box

RETAINER BOXES/MOUTH GUARD BOXES/DENTURE BATHS




600 Retainer Box W75 mm x D61 mm x H28 mm
601 Mouth Guard Box W77 mm x D64 mm x H40 mm
602 All Guard Box W77 mm x D77 mm x H38 mm
603 Denture Bath W83 mm x D68 mm x H52 mm

Colours: Blue, Green, Yellow, White, Lime, Light Blue, Lavender or Assorted

600 Retainer Boxes, 10 pcs/box
601 Mouth Guard Boxes, 10 pcs/box
602 All Guard Boxes, 10 pcs/box
603 Denture Baths, 10 pcs/box

DENTAL APPLIANCE BOXES WITH MIRRORS




PASTEL COLOURS	
Vanilla	Sea Green
Sand	Lilac
Light Mauve	Baby Blue
White	Assorted

- Sturdy strong plastic
- Size: W80 mm x D72 mm x H32 mm

612
Dental Appliance Boxes
10 pcs/box

INTRAORAL CHEEK RETRACTORS



- Holds mouth open, depresses tongue movement
- Anatomically designed
- Available in both autoclavable and disposable



Autoclavable Intraoral Cheek Retractor
2060A Adult Large, 1 pc/box
2060B Adult Small, 1 pc/box
2060C Child, 1 pc/box

Disposable Intraoral Cheek Retractors
2080A Adult Large, 10 pcs/box
2080B Adult Small, 10 pcs/box
2080C Child, 10 pcs/box

INTRAORAL MOUTH LIPPERS



- Expands the mouth vertically rather than horizontally
- Ideal to working in the gingival areas for anteriors, ortho procedures, bleaching etc...
- Available in both autoclavable and disposable



Autoclavable Intraoral Lippers, 2 pcs/box
2040A Adult, Semi-Transparent
2040C Child, Semi-Transparent


Disposable Intraoral Lippers, 10 pcs/box
2030A Adult, White
2030C Child, White

AUTOCLAVABLE CHEEK RETRACTORS



- Autoclavable up to 134°C

Mouth Retractor
2050C Child, 2 pcs/pack
2050A Adult, 2 pcs/pack
2055C Child, 2 pcs/pack
2055A Adult, 2 pcs/pack



AUTOCLAVABLE MOUTH PROPS LATEX FREE

2090A Mouth Prop Adult
2090P Mouth Prop Pado

- Made of chemical resistant silicone
- Autoclavable up to 134°C
- Available in: - Adult (30 mm x 40 mm x 20 mm)
Pado (30 mm x 25 mm x 20 mm)



DISPOSABLE MULTI-PURPOSE DAPPEN DISHES

004B Standard Pack
Assorted Colours
1000 pcs/box

- Disposable to save time and prevent cross-contamination
- For prophylaxis, etchant and bonding etc
- Colours: Lavender, White, Green, Light blue and Yellow

DISPOSABLE PROPHY RINGS



- Disposable to save time and prevent cross-contamination
- For prophylaxis, endo etc
- Colour: Blue

010 Disposable Prophylaxis Rings, 200 pcs/box

AUTOCLAVABLE PROPHY RINGS



- Autoclavable up to 134°C
- For prophylaxis, endo etc
- Colour: Grey

A010 Autoclavable Prophylaxis Rings, 10 pcs/pack

GLASS DAPPEN DISH



- Colours: Light Blue, Green or Clear

040 Glass Dappen Dish

DOUBLE ARCH



LARGE blue



MEDIUM yellow



SMALL green

SINGLE ARCH



LARGE blue



MEDIUM yellow



SMALL green


DISPOSABLE FLUORIDE TRAYS

- Extremely comfortable, anatomically designed trays
- Economical price
- Colour coded

Double Arch, 50 pcs/box
7800S Small (green)
7800M Medium (yellow)
7800L Large (blue)

Single Arch, 100 pcs/box
7801S Small (green)
7801M Medium (yellow)
7801L Large (blue)

DISPOSABLE MIXING PADS



3220

3225

Mixing Pads, 3X100 pcs/pack
3220 Small 35 mm X 45 mm
3225 Large 70 mm X 80 mm

DISPOSABLE MIXING WELLS



906

907

- Sturdy design

906 Mixing Wells 4 Wells 200 pcs/box
907 Mixing Wells 2 Wells 500 pcs/box

actual size



081 ETCH-TIP



080 FLOW-TIP



082 LINER-TIP

DISPOSABLE FLOWABLE COMPOSITE AND ETCHING GEL TIPS

- Universal Luer-lock design
- Bendable
- Gauges: 080 G.20 (0.9 mm)
081 G.23 (0.6 mm)
082 G.18 (1.2 mm)

100 pcs/pack
080 Flow-Tips Black G.20
081 Etch-Tips Light Blue G.23
082 Liner-Tips Pink G.18


DISPOSABLE VALUSTICK STICKY APPLICATORS



- Bends into any angle for placement of Veneers, Inlays, Onlays, Crowns etc
- Strong adhesive for multiple applications
- Disposable
- Length: 110 mm
- Colour: White

909 Valustick Sticky Applicators
50 pcs/box

DISPOSABLE PLASTIC WEDGES



7650X (17 mm)

7650M (13 mm)

7650L (15 mm)

7650S (11 mm)

- Anatomically designed
- Soft tough plastic

Disposable Plastic Wedges, 100 pcs/box
7650X X-Large (17 mm) White
7650L Large (15 mm) Yellow
7650M Medium (13 mm) Green
7650S Small (11 mm) Blue



TONGUE CLEANERS

- Effective in the treatment of halitosis
- Multiple cleaning ridges
- Assorted colours

9000 Tongue Cleaners
20 pcs/pack

OBLONG ARTICULATING MOUNTING PLATES



- Fits Denar type articulators
 - Fully disposable
 - Colours: Blue or White
- 490
Mounting Plates
100 pcs/box


AUTOCLAVABLE ACRYLIC MIXING CUPS



- Autoclavable up to 134°C


280 Mixing Cup Small dia. 25 mm
281 Mixing Cup Medium dia. 35 mm
282 Mixing Cup Large dia. 50 mm
283 Mixing Cup X-Large dia. 70 mm

DISPOSABLE EVACUATION TRAPS



041 Disposable Evacuator Traps dia. 48 mm Green, 20 pcs/pack
042 Disposable Evacuator Traps dia. 54 mm Blue, 20 pcs/pack

AUTOCLAVABLE ALUMINIUM BUR BLOCKS




- Ideal for storage of dental burs
- Specially designed for ultrasonic cleaning and autoclaving of burs
- Sturdy construction, autoclavable up to 134°C
- Colours: Silver and Blue

Small Bur Block
038-FG 15 FG W58 mm x D25 mm x H31 mm
039-MIX 10 FG + 5 RA W58 mm x D25 mm x H31 mm
037-HP 6 HP W54 mm x D28 mm x H62 mm

Large Bur Block
2832-FG 32 FG W80 mm x D42 mm x H33 mm
2832-MIX 16 FG + 16 RA W80 mm x D42 mm x H33 mm
2880-FG 80 FG W100 mm x D75 mm x H33 mm
2880-MIX 40 FG + 40 RA W100 mm x D75 mm x H33 mm


DISPOSABLE EVACUATION TRAPS



- Fits Adec type unit

044 Evacuation Traps
20 pcs/pack

AUTOCLAVABLE PLASTIC BUR DISPENSERS



- Specially designed for cleaning and autoclaving of burs
- Locks open or closed for processing
- Autoclavable up to 134°C
- Size: 75 mm x 35 mm x 10 mm

050-FG, Bur Dispenser 12 FG
050-MIX Bur Dispenser 6 FG + 6 RA

COLOURS	Yellow
Blue	White
Light Blue	Light Beige
Green	Lavender

AUTOCLAVABLE MAGNETIC BUR STANDS



COLOURS	Yellow
Blue	White
Light Blue	Light Beige
Green	Lavender

- Stands autoclavable up to 134°C
- Fits both FG and RA burs
- Size: Small Square 42 mm x 42 mm x 12 mm
Large Rectangular 80 mm x 42 mm x 12 mm

033 Bur Stand Square, 8 Holes 034 Bur Stand Rectangular, 14 Holes
033A Bur Stand Cover Square 034A Bur Stand Cover Rectangular

AUTOCLAVABLE BRASS BUR BRUSH



- Telescopic handle with brass wire bristles
- Length: 12 cm

1110
Brass Bur Brush

DISPOSABLE CURVED TIP SYRINGES



078
Curved Tip
Syringes
50 pcs/box

- Made of polypropylene plastic
- Full 12 ml capacity
- Self sealing plunger
- Designed for perio irrigation, impression materials, post surgical home use, fluoride treatments etc

BRASS / STAINLESS STEEL WIRE BRUSHES

- Autoclavable up to 134°C
- Brass / stainless steel wire brushes remove extra-stubborn stains from the tightest of spaces
- Non-slip handle



1115 Brass Wire Brush
1118 Stainless Steel Wire Brush

AUTOCLAVABLE FRACTURE DETECTOR


The instrument for locating cracks in teeth



- Autoclavable to 134°C
- Simple tool for practice
- Quick and precise to locate cracks
- Can be used for crown and bridge seating


025 Autoclavable
Fracture Detector Blue

GLASS CEMENT MIXING SLAB



- Polished glass mixing slab with bevelled edges
- Size: 65 mm x 135 mm x 15 mm

045 Glass Slab



DISPOSABLE SELF-SEALING STERILISATION POUCHES

- With internal and external steam indicators
- Triple adhesive seal to ensure sterile storage
- Blue tinted plastic window helps detect puncture holes that may be caused by jagged instruments

Standard Pack, 200 pcs/box

- 6001 2 1/4" x 5" (57 mm x 130 mm)
- 6003 3 1/2" x 5 3/4" (90 mm x 145 mm)
- 6005 2 3/4" x 9" (70 mm x 230 mm)
- 6007 3 1/2" x 9" (90 mm x 230 mm)
- 6009 5 1/4" x 10" (135 mm x 260 mm)
- 6011 7 1/2" x 13" (190 mm x 330 mm)

Standard Pack, 100 pcs/box

- 6015 10" x 15" (255 mm x 380 mm)


DISPOSABLE IIR FACE MASKS FLUID RESISTANT LATEX FREE



PERFORMANCE TYPE IIR
EN14683
BFE ≥99% PFE ≥99%
LEVEL 3 BARRIER
ASTM F2100

- Non-irritating
- Fluid resistant, excellent filtration (BFE 99%)
- Permits air to pass through mask to prevent glasses from fogging
- Colours: Blue, Green, Pink, Lavender or Black (ear-loop only)

IIR Face Masks, 50 pcs/box
2288 Ear-Loop
2288T Tie-On



MASK EXTENDERS

- Relieve ear fatigue and pain
- Adjustable buckle strip for earloop masks
- Autoclavable soft silicone
- Colours: Assorted

098
Mask Extenders
5 pcs/pack

DISPOSABLE IIR FACE MASKS WITH SHIELDS EAR-LOOP FLUID RESISTANT LATEX FREE


PERFORMANCE TYPE IIR
EN14683
BFE ≥99% PFE ≥99%
LEVEL 3 BARRIER
ASTM F2100



- Non-irritating
- Fluid resistant, excellent filtration (BFE 99%)
- Anti-fog shields
- Permits air to pass through mask to prevent shield from fogging
- Colours: Blue, Green, Pink or Lavender

FS2288 IIR Face Masks with Shields
Ear-Loop, 25 pcs/box

X-RAY BITE CUSHIONS



- Peel and stick
- Comfortable polyethylene foam
- For both PSP envelopes or sensors with sleeves
- Latex free



HORIZONTAL VERTICAL PERIAPICAL

2018 X-Ray Bite Cushions, 200 pcs/box


DISPOSABLE EYE SHIELDS



- Available in tinted or clear lenses
- Light weight, optical plastic shields
- Assorted colours frames: Blue, Light Blue, Lavender, Green and White

Disposable Eye Shield
20 Lenses + 20 Frames (Assorted Colours)
720T, Tinted
720C, Clear


Replacement Eye Shield Lens Only, 50 pcs/pack
720T-R, Tinted
720C-R, Clear



DISPOSABLE BITE WING LOOPS

- Convenient pop-open design ensures quick and easy use
- Manufactured from quality ten point paper stock
- Quick, simple and economical


Bite Wing Loops, 500 pcs/box
A. 2011-1 Film Child
B. 2011-2 Film Adult
C. 2012-1 Sensor Size 1
D. 2012-2 Sensor Size 2



AUTOCLAVABLE PLASTIC X-RAY FILM HOLDER

- One end for anterior x-rays, the other end for posterior
- Both ends securely hold film without bending it
- Improves quality and accuracy of intraoral films
- Helps to standardize the bisecting angle technique
- Colours: Blue and Light Beige

1188 X-Ray Film Holder



DELUXE ACRYLIC COMPOSITE ORGANISERS

- Smart, flexible design
- Sturdy, durable construction
- Clear material makes contents easy to see
- For storage of composites, adhesive, micro applicators etc.

Composite Organiser
 M060 Mini *W19.5 cm x H9.5 cm x D14.5 cm*
 S060 Standard *W29 cm x H9.5 cm x D14.5 cm*


COTTON ROLL DISPENSER



- Top flips back for easy loading, sliding drawer dispenses up to 4 rolls with a single pull
- Foam rubber strips inserted into base prevent slipping on table top or counter

1140 Roll Dispenser

COLOURS	Blue	Light Blue	Green
Yellow	White	Light Beige	Lavender




DISPOSABLE MICRO APPLICATOR DISPENSERS

Applicator Dispenser
 18-900 Standard Handle
 18-9900 Mini Handle

- Small, sturdy, durable construction

ACRYLIC ORTHO CHAIN DISPENSER



- Long lasting dispenser that will fit the decor of any office
- Chain is dispensed easily with an interior roller
- Holds up to 15 spools


059 Chain Dispenser



ACRYLIC COMPOSITE DISPENSER

- Smart, durable construction
- Clear material makes contents easy to see
- For storage of 20 composite syringes, etch and bonding material etc.

H060 Composite Dispenser



ACRYLIC COMPOSITE SYRINGE HOLDERS

551 Syringe Holder *Small*
 552 Syringe Holder *Medium*

- Smart, durable construction
- Size: 551 Small (11 holes) *W10 cm x D12.5 cm x H12 cm*
 552 Medium (23 holes) *W20 cm x D12.5 cm x H12 cm*



LIGHT PROTECTION HAND SHIELD

- Hand held shield for curing light
- Cold disinfection only
- Size: *12.5 cm x 26 cm x 0.3 cm*

063 Light Protection Hand Shield



LIGHT PROTECTION COVERS

- Extends working time of composite materials
- Sturdy, durable construction
- Cold disinfection only
- Size: Small *W10 cm x D7.5 cm x H3 cm*
 Medium *W15 cm x D10 cm x H3 cm*

062-0 Light Protection Cover *Small*
 062-1 Light Protection Cover *Medium*



ACRYLIC IMPRESSION GUN HOLDERS

- For both loaded or empty guns
 - Store all guns in upright position
 - Include storage compartments
 - Sturdy, durable construction
 - Sizes:
 - 072 W145 mm x D140 mm x H150 mm
 - 074 W265 mm x D140 mm x H150 mm
 - 076 W395 mm x D140 mm x H150 mm
- 072 Gun Holder for 2 guns
074 Gun Holder for 4 guns
076 Gun Holder for 6 guns



ACRYLIC GRINDING BOX

- Ideal for grinding and sand blasting
 - Retains flying debris
 - Sturdy durable construction
 - Size: W30 cm x D20 cm x H18 cm
- 065 Grinding Box

ACRYLIC GLOVE AND TISSUE DISPENSERS



- Smart and attractive acrylic design
 - Sturdy, durable construction
 - Pre-drilled holes for wall mounting
 - Internal measurements:-
 - 064 Glove Dispenser **Standard** W25 cm x D14 cm x H9 cm
 - 066 Glove Dispenser **Small** W25 cm x D14 cm x H7 cm
 - 068 Glove/Tissue Dispenser **Large** W24 cm x D12.5 cm x H10 cm
 - 069 Glove/Tissue Dispenser **Extra Large** W24 cm x D12.5 cm x H11 cm
- 064 Glove Dispenser **Standard**
066 Glove Dispenser **Small**
068 Glove/Tissue Dispenser **Large**
069 Glove/Tissue Dispenser **Extra Large**




PLASTIC CUP DISPENSER

- Holds approx. 50 disposable cups from 150 ml to 200 ml (5 oz to 7 oz)
- Self adhesive
- Made of ABS plastic (acrylonitrile butadiene styrene)
- Size: dia. 100 mm x L280 mm

1280 Plastic Cup Dispenser

ACRYLIC BIB AND C-FOLD TOWEL DISPENSERS




- Sturdy, durable construction
 - Smart, attractive acrylic design
 - Sizes:
 - 057 W15 cm x D11 cm x H12 cm
 - 058 W23.5 cm x D10 cm x H13.5 cm
- 057 Bib Dispenser
058 C-Fold Towel Dispenser
057CV Bib Dispenser With Cover
058CV C-Fold Towel Dispenser With Cover



ACRYLIC MASK DISPENSER

095 Mask Dispenser

- Sturdy, durable construction
- Hinged lid for protection against debris
- Clear material makes contents easy to see
- Pre-drilled holes for wall mounting
- Size: W194 mm x D106 mm x H108 mm




ACRYLIC POUCH DISPENSERS

- Can be used on table top or wall mounted in either vertical or horizontal position
- Sizes:-
 - 061-0 (7.5 cm x 14.5 cm x 21 cm)
 - 061-1 (8.5 cm x 27.5 cm x 21 cm)
 - 061-2 (10 cm x 18.5 cm x 21 cm)
 - 061-3 (10 cm x 27.5 cm x 21 cm)
 - 061-4 (15 cm x 29.5 cm x 21 cm)

061-0 Pouch Dispenser for 5.7 cm x 10 cm Pouches
061-1 Pouch Dispenser for 7 cm x 23 cm Pouches
061-2 Pouch Dispenser for 9 cm x 13.5 cm Pouches
061-3 Pouch Dispenser for 9 cm x 23 cm Pouches
061-4 Pouch Dispenser for 13.5 cm x 26 cm Pouches


ACRYLIC BARRIER FILM DISPENSER



- Sturdy durable construction
- Clear material makes contents easy to see
- Size: dia. W14 cm x D23 cm x H19 cm

1155 Barrier Film Dispenser

CROSS-SECTION



E01 SUPER FILES GOLD NiTi heat activation

- Smart varied taper design
- High cutting efficiency
- Good flexibility and fracture resistance
- Torque 2.0 Ncm - 2.5 Ncm
- Recommended speed 150-300 r.p.m.


E01 SUPER FILES GOLD NiTi heat activation			
Length	Size	Taper %	Tip (mm)
19 mm	SX	3.5-9	0.19
21 mm or 25 mm	S1	2-11	0.17
	S2	4-11.5	0.20
	F1	7	0.20
	F2	8	0.25
	F3	9	0.30
	F4	6	0.40
	F5	5	0.50

E01 Super Files Gold
Single or Assorted Sizes
6 pcs/pack

E01 RECOMMENDED OPERATION PROCESS

PROBE	#10	
	#15	
OPENING	SX	
SHAPING	S1	
	S2	
FINISHING	F1	
	F2	
	F3	
ENLARGING	F4	
	F5	

CROSS-SECTION



E02 SUPER FILES BLUE NiTi heat activation

- Smart varied taper design
- High flexibility and pre-curved memory
- Safe tip, reduces the possibility of root canal deviation
- Torque 2.0 Ncm - 2.5 Ncm
- Recommended speed 150-300 r.p.m.


E02 SUPER FILES BLUE NiTi heat activation			
Length	Size	Taper %	Tip (mm)
19 mm	SX	3.5-9	0.19
21 mm or 25 mm	S1	2-11	0.17
	S2	4-11.5	0.20
	F1	7	0.20
	F2	8	0.25
	F3	9	0.30
	F4	6	0.40
	F5	5	0.50

E02 Super Files Blue
Single or Assorted Sizes
6 pcs/pack

E02 RECOMMENDED OPERATION PROCESS

PROBE	#10	
	#15	
OPENING	SX	
SHAPING	S1	
	S2	
FINISHING	F1	
	F2	
	F3	
ENLARGING	F4	
	F5	

CROSS-SECTION



- Simple effective system
- Unique flexibility
- Torque 2.0 Ncm - 2.5 Ncm
- Recommended speed 150-300 r.p.m.

E03 Pro-Flexi Files
Single or Assorted Sizes
6 pcs/pack


E03 PRO-FLEXI FILES NiTi heat activation function

E03 PRO-FLEXI FILES heat activation function			
Length	Size	Taper %	Tip (mm)
19 mm	#01/09	9	0.25
21 mm or 25 mm	#20/04	4	0.20
	#25/04	4	0.25
	#30/04	4	0.30
	#35/04	4	0.35
	#25/06	6	0.25

E03 RECOMMENDED OPERATION PROCESS

OPENING	#01	
NARROW CANAL	#20 (-4%)	
	#25 (-4%)	
MEDIUM CANAL	#30 (-4%)	
WIDE CANAL	#35 (-4%)	
	#25 (-6%)	

CROSS-SECTION



E05 ONE FILES NiTi heat activation

- Reciprocating movement
- Use only one instrument for root canal preparation, efficient and time saving
- 3 sizes for different size root canals
- Torque 2.0 Ncm - 2.5 Ncm
- Recommended speed 150-300 r.p.m.

E05 One Files
Single or Assorted Sizes
6 pcs/pack

E05 ONE FILES NiTi heat activation	Length	Size	Taper %	Tip (mm)
	21 mm or 25 mm	#R25	8	0.25
		#R40	6	0.40
		#R50	5	0.50

E05 ONE FILES	#R25	
	#R40	
	#R50	

E38 BARBED BROACHES stainless steel



- Good flexibility
- Highly resistant to fracture

E38 Barbed Broaches
Single Size
Stainless Steel
10 pcs/pack

E38 BARBED BROACHES stainless steel	
Length	Size
21 mm	XXXXF (0)
	XXXF (1)
	XXF (2)
	XF (3)
	F (4)
	M (5)
	C (6)

E48 C+ FILES stainless steel

- For calcified root canal, with 3 mm very rigid tip
- Recommended for single use only

E48 C+ Files
Single or Assorted Sizes
6 pcs/pack

E48 C+ FILES stainless steel	
Length	Size
21 mm or 25 mm	#06
	#08
	#10
	#12.5
	#15

E04 ULTRASONIC FILES NiTi alloy

- For use with ultrasonic scaler
- Excellent tenacity
- Length: 33 mm



E04 ULTRASONIC FILES NiTi alloy	
#15	#30
#20	#35
#25	#40

E04
Ultrasonic Files
Single or Assorted Sizes
6 pcs/pack

E11 SPIRAL FILLERS stainless steel

E11 SPIRAL FILLERS stainless steel	
Length	Size
21 mm or 25 mm	#25
	#30
	#35
	#40

E11
Spiral Fillers
Single or Assorted Sizes
4 pcs/pack

E12 GATES DRILLS stainless steel

- High cutting efficiency
- Recommended speed 800-1,200 r.p.m.



E12 GATES DRILLS stainless steel		
Length	Size	Diameter (mm)
28 mm or 32 mm	#1	0.5
	#2	0.7
	#3	0.9
	#4	1.1
	#5	1.3
	#6	1.5

E12
Gates Drills
Single or Assorted Sizes
6 pcs/pack

E15 PEESO REAMERS stainless steel

- High cutting efficiency
- Recommended speed 800-1,200 r.p.m.



E15 PEESO REAMERS stainless steel		
Length	Size	Diameter (mm)
28 mm or 32 mm	#1	0.7
	#2	0.9
	#3	1.1
	#4	1.3
	#5	1.5
	#6	1.7

E15
Peeso Reamers
Single or Assorted Sizes
6 pcs/pack

E30 K-FILES stainless steel

E30 K-FILES stainless steel					
#06	#08	#10			
#15	#20	#25	#30	#35	#40
#45	#50	#55	#60	#70	#80

E30 K-Files Stainless Steel 21 or 25 mm, 6 pcs/pack
No.06-80 Single Size
No.15-40, 45-80 Assorted Sizes

E32 K-FILES NiTi alloy

E32 K-FILES NiTi alloy					
#08	#10				
#15	#20	#25	#30	#35	#40
#45	#50	#55	#60	#70	#80

E32 K-Files NiTi Alloy 21 or 25 mm, 6 pcs/pack
No.06-80 Single Size
No.15-40, 45-80 Assorted Sizes

E34 H-FILES stainless steel

E34 H-FILES stainless steel					
#08	#10				
#15	#20	#25	#30	#35	#40
#45	#50	#55	#60	#70	#80

E34 H-Files Stainless Steel 21 or 25 mm, 6 pcs/pack
No.10-80 Single Size
No.15-40, 45-80 Assorted Sizes

E36 H-FILES NiTi alloy

E36 H-FILES NiTi alloy					
#08	#10				
#15	#20	#25	#30	#35	#40
#45	#50	#55	#60	#70	#80

E36 H-Files NiTi Alloy 21 or 25 mm, 6 pcs/pack
No.10-80 Single Size
No.15-40, 45-80 Assorted Sizes

E40 / E42 SPREADERS stainless steel / NiTi alloy

E40 SPREADERS stainless steel					
#15	#20	#25	#30	#35	#40

E40 Spreaders Stainless Steel 21 or 25 mm, 6 pcs/pack
No.15-40, Single or Assorted Sizes

E42 SPREADERS NiTi alloy					
#15	#20	#25	#30	#35	#40

E42 Spreaders NiTi Alloy 21 or 25 mm, 6 pcs/pack
No.15-40, Single or Assorted Sizes

E44 / E46 PLUGGERS stainless steel / NiTi alloy

E44 PLUGGERS stainless steel					
#15	#20	#25	#30	#35	#40

E44 Pluggers Stainless Steel 21 or 25 mm, 6 pcs/pack
No.15-40, Single or Assorted Sizes

E46 PLUGGERS NiTi alloy					
#15	#20	#25	#30	#35	#40

E46 Pluggers NiTi Alloy 21 or 25 mm, 6 pcs/pack
No.15-40, Single or Assorted Sizes




REFINE MaxPiezo 7+ / DS7+ PORTABLE LED SCALERS

- Automatic frequency tracking, real-time feedback
- Extra large, capacity 600 ml water tank for water, chlorhexidine, hydrogen peroxide, sodium hypochlorite etc...
- Power and water controlled by touch panel, easy to clean for infection control
- Detachable handpiece
- Frequency: **28 kHz ± 3 kHz**
- Complete with 8 tips
- Size: **W180 mm x D180 mm x H150 mm**
- Weight: **1.6 Kg (unit only)**
- Colour: White

ONE YEAR WARRANTY

(except scaler tips)

MaxPiezo 7 (EMS Type) Portable **Non-LED** Scaler Set
 MaxPiezo 7+ (EMS Type) Portable **LED** Scaler Set
 MaxPiezo DS7 (SATELEC Type) Portable **Non-LED** Scaler Set
 MaxPiezo DS7+ (SATELEC Type) Portable **LED** Scaler Set




REFINE MaxPiezo 6+ / DS6+ LED SCALERS

- Automatic frequency tracking, real-time feedback
- Power and water controlled by touch panel, easy to clean for infection control
- Detachable handpiece
- Frequency: **28 kHz ± 3 kHz**
- Complete with 6 tips
- Size: **W155 mm x D168 mm x H80 mm**
- Weight: **0.6 Kg (unit only)**
- Colour: White

ONE YEAR WARRANTY

(except scaler tips)

MaxPiezo 6 (EMS Type) **Non-LED** Scaler Set
 MaxPiezo 6+ (EMS Type) **LED** Scaler Set
 MaxPiezo DS6 (SATELEC Type) **Non-LED** Scaler Set
 MaxPiezo DS6+ (SATELEC Type) **LED** Scaler Set




REFINE N3000L / V3000L BUILT-IN LED SCALERS

- Automatic frequency tracking, real-time feedback
- Detachable handpiece
- Frequency: **28 KHz ± 3 KHz**
- Complete with 6 tips

ONE YEAR WARRANTY

(except scaler tips)

N3000 (EMS Type) Built-In **Non-LED** Scaler Set
 N3000L (EMS Type) Built-In **LED** Scaler Set
 V3000 (SATELEC Type) Built-In **Non-LED** Scaler Set
 V3000L (SATELEC Type) Built-In **LED** Scaler Set



REFINE AUTOCLAVABLE SCALER HANDPIECE KITS

- Complete with:-
 1x Handpiece 3H / 5L / 7H / 7L
 3x Perio scaler tip P1 / PD1
 1x Torque wrench TW-1L
 1x Scaler tip holder No.1
 1x Sterilization box



3H **Non-LED** Handpiece Kit (EMS Type)
 7H **Non-LED** Handpiece Kit (SATELEC Type)
 5L **LED** Handpiece Kit (EMS Type)
 7L **LED** Handpiece Kit (SATELEC Type)
 (for handpiece connections, please refer to Page 37)

ONE YEAR WARRANTY

(except scaler tips)

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France

Scaler Handpieces and Accessories

REFINE SCALER HANDPIECES



HP-3H Non-LED Scaler Handpiece (EMS Type)



HS-7H Non-LED Scaler Handpiece (SATELEC Type)



HP-5L LED Scaler Handpiece (EMS Type)



HS-7L LED Scaler Handpiece (SATELEC Type)

REFINE SCALER TORQUE WRENCHES



TW-1L
For EMS¹,
KAVO PiezoLED⁵
SATELEC², NSK³
type scalers



TS-1L
For SIRONA⁴
type scaler



TK-1L
For KAVO⁵
SONICflex⁵
SONOsoft⁵
type scalers

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
2. SATELEC is a registered trademark of Satelec/Acteon Group France
3. NSK is a registered trademark of Nakanishi Inc. Japan
4. SIRONA is a registered trademark of Siroa Dental System
5. KAVO, SONICflex, SONOsoft, PiezoLED are registered trademarks of Kavo Dental GmbH Germany

REFINE SCALER HANDPIECE HOSES




3H Non-LED (EMS Type)
Handpiece Hose



7H Non-LED (SATELEC Type)
Handpiece Hose



5L LED (EMS Type)
Handpiece Hose



7L LED (SATELEC Type)
Handpiece Hose



E04 ULTRASONIC FILES NiTi alloy

- For use with ultrasonic scaler
- Excellent tenacity
- Length: 33 mm

E04 ULTRASONIC FILES NiTi alloy	
#15	#30
#20	#35
#25	#40


E04 Ultrasonic Files
Single or
Assorted Sizes
6 pcs/pack

REFINE SCALER TIP HOLDER



NO.2

REFINE ENDO WRENCH



Multi function wrench
for both endo chuck & bur

TW-A

REFINE FOOT SWITCH




REFINE SCALER TIPS STAINLESS STEEL

FOR EMS¹, SATELEC², NSK³ and SIRONA⁴ REFINE Type Scalers


SCALING




G1 (EMS¹ Type)
GD1 (SATELEC², NSK³ Type)
GS1 (SIRONA⁴ Type)
 Supragingival deposits in all quadrants




G2 (EMS¹ Type)
GD2 (SATELEC², NSK³ Type)
GS2 (SIRONA⁴ Type)
 Heavy supragingival deposits on anterior teeth




G3 (EMS¹ Type)
GD3 (SATELEC², NSK³ Type)
GS3 (SIRONA⁴ Type)
 Supragingival deposits in all quadrants, including the interproximal and sulcus areas



G4 (EMS¹ Type)
GD4 (SATELEC², NSK³ Type)
GS4 (SIRONA⁴ Type)
 Supragingival and interproximal deposits




G5 (EMS¹ Type)
GD5 (SATELEC², NSK³ Type)
GS5 (SIRONA⁴ Type)
 Supragingival deposits in all quadrants




G6 (EMS¹ Type)
GD6 (SATELEC², NSK³ Type)
GS6 (SIRONA⁴ Type)
 Voluminous supragingival deposits


ENDO



E1 (EMS¹ Type)
ED1 (SATELEC², NSK³ Type)
ES1 (SIRONA⁴ Type)
 120°file holder for anterior root canal irrigation




E2 (EMS¹ Type)
ED2 (SATELEC², NSK³ Type)
ES2 (SIRONA⁴ Type)
 95°file holder for molar teeth root canal irrigation




E3 (EMS¹ Type)
ED3 (SATELEC², NSK³ Type)
ES3 (SIRONA⁴ Type)
 Gutta-percha melting




E3D (EMS¹ Type)
ED3D (SATELEC², NSK³ Type)
ES3D (SIRONA⁴ Type)
 Diamond coated instrument for root canal enlarging




P4D (EMS¹ Type)
PD4D (SATELEC², NSK³ Type)
PS4D (SIRONA⁴ Type)
 Diamond coated instrument for root canal enlarging




E14 (EMS¹ Type)
ED14 (SATELEC², NSK³ Type)
ES14 (SIRONA⁴ Type)
 Removing root filling obstructions during root canal retreatments




E15 (EMS¹ Type)
ED15 (SATELEC², NSK³ Type)
ES15 (SIRONA⁴ Type)
 Removing root filling obstructions during root canal retreatments



E4 (EMS¹ Type)
ED4 (SATELEC², NSK³ Type)
ES4 (SIRONA⁴ Type)
 Removing root filling obstructions during root canal retreatment



E4D (EMS¹ Type)
ED4D (SATELEC², NSK³ Type)
ES4D (SIRONA⁴ Type)
 Removing calcification and filling material during root canal retreatment



E5 (EMS¹ Type)
ED5 (SATELEC², NSK³ Type)
ES5 (SIRONA⁴ Type)
 Removing root filling obstructions during root canal retreatment



E5D (EMS¹ Type)
ED5D (SATELEC², NSK³ Type)
ES5D (SIRONA⁴ Type)
 Removing calcification and filling material during canal retreatment


1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France
 3. NSK is a registered trademark of Nakanishi Inc. Japan
 4. SIRONA is a registered trademark of Sirona Dental System

REFINE SCALER TIPS STAINLESS STEEL


FOR EMS¹, SATELEC², NSK³ and SIRONA⁴ Type Scalers




P1 (EMS¹ Type)
PD1 (SATELEC², NSK³ Type)
PS1 (SIRONA⁴ Type)
 Shallow pocket subgingival deposits



P3 (EMS¹ Type)
PD3 (SATELEC², NSK³ Type)
PS3 (SIRONA⁴ Type)
 Deep pocket subgingival deposits



P4 (EMS¹ Type)
PD4 (SATELEC², NSK³ Type)
PS4 (SIRONA⁴ Type)
 Shallow pocket subgingival deposits



P3D (EMS¹ Type)
PD3D (SATELEC², NSK³ Type)
PS3D (SIRONA⁴ Type)
 Diamond coated instrument for root planning after periodontal flap surgery

PERIO




P2L (EMS¹ Type)
PD2L (SATELEC², NSK³ Type)
 Left angle instrument to remove calculus from very narrow inter-root space and furcation



P2R (EMS¹ Type)
PD2R (SATELEC², NSK³ Type)
 Right angle instrument to remove calculus from very narrow inter-root space and furcation



P2LD (EMS¹ Type)
PD2LD (SATELEC², NSK³ Type)
 Removal of calculus from narrow inter-root spaces (diamond-coated)



P2RD (EMS¹ Type)
PD2RD (SATELEC², NSK³ Type)
 Diamond coated instrument for root planning

FOR AMDENT⁶ Type Scalers



A1 (AMDENT⁶ Type)
 Supragingival scaling (#37)



A2 (AMDENT⁶ Type)
 Supragingival scaling (#39)

FOR KAVO SONICflex⁵ Type Scalers



GK1 (KAVO⁵ Type)
 Supragingival deposits in all quadrants



GK2 (KAVO⁵ Type)
 Supragingival heavy calculus and plaque




GK3 (KAVO⁵ Type)
 Supragingival interdental calculus and calculus at the neck of the tooth




GK4 (KAVO⁵ Type)
 Supragingival deposits, including the interproximal and sulcus areas



GK5 (KAVO⁵ Type)
 Subgingival calculus




GK6 (KAVO⁵ Type)
 Left angle instrument for removing subgingival calculus




GK7 (KAVO⁵ Type)
 Right angle instrument for removing subgingival calculus

FOR KAVO SONOsoft⁵ Type Scalers



GC1 (KAVO⁵ Type)
 Supragingival deposits in all quadrants




GC2 (KAVO⁵ Type)
 Supragingival interdental calculus and calculus at the neck of the tooth



PC1 (KAVO⁵ Type)
 Subgingival calculus, including the interproximal and sulcus areas


KAVO PiezoLED⁵ Type Scalers



201 (KAVO⁵ Type)
 Supragingival and interproximal deposits



202 (KAVO⁵ Type)
 Supragingival deposits in all quadrants, including the interproximal and sulcus areas



203 (KAVO⁵ Type)
 Deep pocket subgingival deposits

1. EMS is a registered trademark of Electro Medical System S.A. Switzerland
 2. SATELEC is a registered trademark of Satelec/Acteon Group France
 3. NSK is a registered trademark of Nakanishi Inc. Japan
 4. SIRONA is a registered trademark of Sirona Dental System
 5. KAVO, SONICflex, SONOsoft, PiezoLED are registered trademarks of Kavo Dental GmbH Germany
 6. AMDENT is a registered trademark of Amdent AB Sweden

Premium Plus International Ltd.
Flat 1001, Yuen Long Trading Centre,
No.33 Wang Yip Street West,
Yuen Long
Hong Kong
+852 2761 3211
mailbox@premiumplushk.com
www.premiumplushk.com

Premium Plus Poland sp. z o.o.
ul. Bukowska 27
62-081
Wysogotowo
Poland
+48 61 880 1094
info@premiumpluspl.com
www.premiumpluspoland.com

Premium Plus UK Ltd.
Unit 2, Knighton Heath Ind. Estate,
Bournemouth,
Dorset BH11 8NE
UK
+44 1202 611011
info@premiumplusuk.com
www.premiumplusuk.com

Premium Plus Dongguan Limited
Flat 101, No.122 Tianqiao Road
Tutang, Changping, Dongguan,
Guandong
China
+86 769 83756768
info@premiumpluschina.com
www.premiumpluschina.com

Premium Plus Dental Supplies Inc.
155 Liberty Avenue
Brooklyn
New York 11212
United States of America
+1 718 366 8686
info@premiumplusdentalusa.com
www.premiumplususa.com

Premium Plus Japan Co. Ltd.
5th Floor, 1-1-31, Nishimidorigaoka,
Toyonakashi
Osaka 560-0005
Japan
+81 6 6845 0066
info@premiumplus.jp
www.premiumplus.jp