

ErgoDenta

Bestseller Catalog

Product Catalog 2025
The Nordic Selection



Comfort in every
instrument


Discover 2500+
Dental Instruments
for Daily Care



At ErgoDenta, we believe strong partnerships lead to success. That's why we work directly with our distributors—no middlemen, no extra costs—just over 2,500 reliable dental instruments delivered straight from our Denmark-based warehouse. From sourcing materials to production and shipping, we handle every step to ensure quality, efficiency, and smooth deliveries.

Inside this catalog, you'll find smart solutions designed to make dental work easier. ErgoRazor scalers and curettes stay sharp up to 3 times longer, reducing replacements. ErgoSlip instruments come in gold, black, and blue coatings, combining durability with a sleek look. ErgoTip allows dentists to replace only the worn-out tip, saving costs and reducing waste. ErgoTrays keep workspaces organized, while ErgoLite offers a cost-effective option without compromising on quality.

With ErgoDenta, you get high-quality instruments, direct from the source. Browse the catalog and find the right products for your market.



DANISH DESIGN

With our Danish roots and skilled craftsmanship, every product is made with great care and attention to detail and quality.

ERGONOMIC FOCUS

We specialize in creating instruments that prioritize comfort and efficiency, helping dental professionals to work better and feel better.

GLOBAL COMPLIANCE

Our products meet the standards of MDR, UKCA, and ISO 13485, ensuring reliable, safe, and top-quality instruments for all.

INNOVATIVE SOLUTIONS

We continuously explore new ways to enhance the dental experience, offering cost-effective tools that support distributors, dentists, and private label needs worldwide.

ErgoRazor Sharp Scalers & Curettes



Upto
**3X Longer
Sharpness**

Long-lasting sharpness and
durability, cryo treated in the UK

Highlights
Read more on Page 44



ErgoSlip

Non-stick composite instruments
Smooth composite work with
black, blue, or gold coatings




Highlights
Read more on Page 30



ErgoTip

interchangeable dental
instrument tips



USE

Perform your procedure with the ErgoTip instrument.


REPLACE

When the tip wears out, simply swap it for a new one.

REPEAT

Keep using the same handle—no need to buy or sharpen the tips anymore.

Highlights
Read more on Page 10



ErgoLite

Ergonomic, Durable & Cost-Effective Instruments



Highlights
Read more on Page 101



ErgoSteel



Ø6, Ø8, Ø9.5 mm options
Anti-rust stainless steel
Secure grip pattern
Light-weight
Ergonomic balance

Handles

ErgoX
Plus

Light
Weight

Perfectly
Balanced

Feel the
Grip

Full Autoclaveable

Anti-Rust
Instrument

Excellent
Craftsmanship

Ergo X



Comes in 15 exciting colors

Handles

TABLE OF CONTENTS

ErgoTip

Page no. 10-11

ErgoRazor

Page no. 44-47

ErgoSlip

Page no. 30-32

ErgoTrays

Page no. 95-97

01

Diagnostics

Explorers	13
Periodontal Probes	14
Explorer & Probes	15
Mouth Mirrors & Handles	17
HD & Endo Mirrors	19
Diagnostic kits	20

02

RESTORATION

Excavators	22
Composite & Filling instruments	24
Composite kits	32
Rubberdam instruments	34
Crown instruments	37
Matrix Retainers	37

03

Peridontics

Gracey Curettes	39
Sickle Scalers	41
Universal Curettes	42
Perio kits	48

04

Extraction

Extraction Forceps	51
Root tip elevators	56

05

Surgery

Scissors	62
Needle holders	68
Hemostatic Forceps	70
Tweezers	71
Retractors	74
Syringes	76
Endo instruments	77
Orthodontics	80
Implantology	84
Bone curettes	87
Periosteal elevators	88
Surgery kits	90

06

Diverse

Sterilization trays	93
Diverse trays & instruments	98
Endo accessories	100
ErgoTape & Color-coded rings	101

07

Ergolite


ErgoLite X instruments	103
ErgoLite S instruments	105
ErgoLite elevators	107

ErgoTip

With ErgoTip, you don't replace the entire instrument—just the worn-out tip. It's that simple.

- **USE:** Perform your procedures with ErgoTip's ergonomic design.
- **REPLACE:** Swap the worn-out tip—no need to replace the entire instrument.
- **REPEAT:** Keep using the same intact handle, saving time, money, and planet resources.


Available in packs of 1 or 5 pieces. The 5-piece pack includes the handle of your choice and includes a wrench




Eco-friendly, cost-effective & efficient

Diagnostics

Explorers




Probes



ErgoTip

With ErgoTip, you don't replace the entire instrument—just the worn-out tip. It's that simple.



Periodontics

Gracey curettes

Gracey 1-2



Gracey 5-6



Gracey 7-8



Gracey 11-12



Gracey 13-14



Mini Gracey 1-2



Mini Gracey 11-12



Mini Gracey 13-14



Sickle Scalers

204S



204SD



M23



Micro Sickle Scaler



Mini Sickle Scaler



H6-H7



ErgoTip handles & accessories

ErgoTip handles available in single-end & double-end. Each handle includes a wrench.



Wrench 5757

Diagnostics




Explorers

Designed for everyday use, these explorers help you check teeth surfaces easily. Available in single and double-end options with comfortable handles.

Single Ended



Double Ended



Explorers

Double Ended



Explorer 17-23

1107



Explorer 17C-23

1108



Weston Explorer

1111



ErgoSteel

ErgoSteel Mini (S)

Length 100mm

Ø6.0mm

ErgoSteel Standard (SS)

Length 104mm

Ø8.0mm

ErgoSteel Plus (SP)

Length 106mm

Ø9.5mm

ErgoX

ErgoX

Length 100mm

Ø11.5mm

ErgoX Plus

Length 100mm

Ø13.6mm

Periodontal Probes

Single Ended

Single and double-end probes with clear markings. Designed for easy and comfortable use.



Probe CP 2
2-4-6-8-10-12mm

1050



Probe CP 2B
Ball End 2-4-6-8-10-12mm

1051



Probe CP 11
3-6-8-11mm

1052



Probe CP 12
3-6-9-12mm

1053



Probe CP 12B
3-6-9-12mm

1054



Probe CP 15
1-15mm

1055



Probe CP 18
3-5-8-10mm

1056



WHO Probe
Ball End 3.5-5.5-8.5-11.5mm

1057



William Probe
3-5-8-11mm

1058



Unlike CP probes that has color coding, William probes are machine-marked

Periodontal Probes

Double Ended



Furcation Probe
1113
without measurements



Nabers Probe
2-4-6-8-10-12mm
1200



**Explorer 6 +
Probe CP 2**
2-4-6-8-10-12mm
1300



**Explorer 6L +
Probe CP 15**
1-15mm
1301



**Probe William +
Explorer Holst**
1112



**Explorer 8 +
WHO Probe**
3.5-5.5-8.5-11.5mm
1304




**Explorer 8 +
Probe CP 15**
1-15mm
1305



**Explorer 9 +
Probe CP 18**
3-5-8-10mm
1306



**Explorer 23 +
Probe CP 2**
2-4-6-8-10-12mm
1307



**Explorer 23 +
Probe CP 2B**
2-4-6-8-10-12mm
1308



**Explorer 23 +
WHO Probe**
3.5-5.5-8.5-11.5mm
1309



ErgoSteel

ErgoSteel Mini (S) Length 100mm Ø6.0mm

ErgoSteel Standard (SS) Length 104mm Ø8.0mm

ErgoSteel Plus (SP) Length 106mm Ø9.5mm

ErgoX

ErgoX Length 100mm Ø11.5mm

ErgoX Plus Length 100mm Ø13.6mm

Periodontal Probes

Double Ended



**Explorer 23 +
Probe CP 12**

3-6-9-12mm

1310



**Explorer 23 +
Probe CP 12B**

3-6-9-12mm

1311



**Explorer Holst +
Probe CP 18**

3-5-8-10mm

1313



**Explorer 23 +
Probe CP 15**

1-15mm

1312




ErgoSteel

ErgoSteel Mini (S) Ø6.0mm



ErgoSteel Standard (SS) Ø8.0mm



ErgoSteel Plus (SP) Ø9.5mm



ErgoX

ErgoX Ø11.5mm



ErgoX Plus Ø13.65mm



A combined dental explorer and probe is very useful during check-ups. The explorer tip finds tiny tooth issues, and the probe measures gum pockets to catch gum disease early. Having both ends in one instrument saves time and streamlines the examination.

Mirror Handles

ErgoSteel Mirror Handles

**ErgoSteel
Small
Mirror handle**



1450 S

Ø 6mm

Length 125mm

**ErgoSteel
Standard
Mirror handle**



1450 SS

Ø 8mm

Length 125mm

**ErgoSteel
Plus
Mirror handle**



1450 SP

Ø 9.5mm

Length 125mm

**Chrome
Plated
Mirror handle**



1453 SS

Ø 8mm

Length 125mm

ErgoX

Mirror handles

Find your perfect fit with our ergonomic mirror handles—15 colors, two sizes

Handles



Mirror, Peach Blush
1401
X | XP



Mirror, Magenta Blue
1402
X | XP



Mirror, Deep Blue
1403
X | XP

15
colors

2
handles



Mirror, Sky Cyan
1404



Mirror, Honey Brown
1405
X | XP



Mirror, Aqua Light
1406
X | XP



Mirror, Petal Pink
1408
X | XP



Mirror, Ruby Pink
1409
X | XP



Mirror, Citrus Orange
1407
X | XP



Mirror, Earth Brown
1411
X | XP



Mirror, Teal Deep
1412
X | XP



Mirror, Lilac Blue
1410
X | XP



Mirror, Lemon Green
1414
X | XP



Mirror, Pure Black
1415
X | XP




Mirror, Misty Cyan
1413
X | XP

Front-Surface

No double image – ensures clear and single reflection.
Suitable for all SS-thread mirror handles

Size	Art nr.
No. 3 - Ø20mm	1503F
No. 4 - Ø22mm	1504F
No. 5 - Ø24mm	1505F

12pcs / box




Plane-Surface

Suitable for all SS-thread mirror handles

Size	Art nr.
No. 3 - Ø20mm	1503P
No. 4 - Ø22mm	1504P
No. 5 - Ø24mm	1505P

12pcs / box



ErgoLite mirror handles

- Light weight, ergonomic mirror handle
- Suitable for Microscope & HD dental mirrors
- Also fits all SS socket dental mirrors
- Available in 10 colors



Colors	Art nr.
Ruby Pink	9309
Lilac Blue	9310
Misty Cyan	9313
Lemon Green	9314
Pure Black	9315
Citrus Orange	9307
Petal Pink	9308
Honey Brown	9305
Sky Cyan	9304
Magenta Blue	9302

HD Mirrors

Experience ultra-clear reflections with our HD mirrors. Their advanced antifog treatment ensures a sharp, single image every time, even during long procedures. Ideal for accurate and reliable dental examinations.

Ultra
clear view

40%
Brighter than FS

Antifog

HD Mirrors in box

Available in 4 colors and 2 sizes
Compatible with ErgoLite and ErgoSteel mirror handles (does not fit with ErgoX or ErgoX Plus mirror handles)



Color	Ø22mm	Ø24mm
Black	1514HK	1515HK
Green	1514HG	1515HG
Blue	1514HB	1515HB
Pink	1514HP	1515HP



Get 2 free ErgoLite mirror handles with every HD mirror box purchase!


Free
handles

HD Microscopic Mirrors

36-layer coated, anti-glare mirror for clear visibility in microscope-assisted procedures. Available in two thicknesses: ultra-thin 0.7mm for maximum access and visibility, and standard 1.5mm for added durability. Sold with the ErgoLite mirror handle.

Regular

Ultrathin



Size	Regular	Ultrathin
3mm	7091R	7091U
5mm	7092R	7092U
7mm	N/A	7095U
2x7mm	7093R	7093U
3x9mm	7094R	7094U



Microscopic mirrors are sold together with ErgoLite mirror handle.



5mm



3x9mm



3mm



2x7mm

Diagnostics Kits


Basic Diagnostic Kit

9908

Available in



X XP SP



Instruments

Dental Mirror Handle
1401

Explorer 5-9
1104

Explorer 23- Periodontal
Probe WHO
1309

Excavator 63-64 (1.5mm)
2802

College Tweezer, Serrated 15cm
2367-1

5 instruments + cassette


Super Diagnostic Kit

9909

Available in



X XP SP



Instruments

Meriam Tweezer, TC
3636-2A

Explorer 23- Periodontal Probe WHO
1107

Explorer 9
1002

Perio Probe CP12
1053

ErgoMirror handle
1401

Plastic Instrument-Posterior
2548

Nystrom III
2525

Excavator 63-64 (1.5mm)
2802

8 instruments + cassette

Restoration



Excavators

Available in round, action, and spoon types with 5 handle options, designed for efficient removal of decayed tissue.

Round Excavators



Round Excavator,
1.0mm

2801



Round Excavator,
1.5mm

2802



Round Excavator,
2.0mm

2803



Round Excavator,
2.5mm

2804



Round Excavator,
3.0mm

2805



Excavator 129/130,
1.7mm

2817



Action Excavators



Action Excavator,
1.0mm

2806



Action Excavator,
1.5mm

2807



Action Excavator,
2.0mm

2808



Action Excavator,
3.5mm

2810



Spoon Excavators



Know your Excavators

Our dental excavators are designed for safe, efficient decay removal

Action Excavator
Combines the benefits of spoon and round types for versatile, effective performance.



Round Excavator
Reaches tight curves for precise cleaning.

Spoon Excavator
Gently scoops out soft decay while preserving healthy tissue.

Composite & Filling instruments

Perfect for shaping, placing, and finishing composites and fillings. Available in various designs and handle options for easy handling and comfort.

Pluggers & Burnishers



Plugger,
0.7mm-1.0mm
2501



Plugger,
1.0mm-2.0mm
2502



Plugger,
2.0mm-2.5mm
2503



Plugger Composite,
1.3mm-1.6mm
2554



Plugger, Cylindrical,
1.2mm-2.0mm
2559



Plugger Clubbed,
1.4mm-1.8mm
2555



Plugger, Cylindrical,
1.8mm-2.5mm
2578



Plugger, Cylindrical,
1.8mm-2.5mm
2561



Plugger Flat
2556



Plugger, Cylindrical,
1.0mm-1.5mm

2560



Pluggers & Burnishers



Plugger, 2.5mm-3.8mm


2581



Burnisher & Plugger,
1.4mm-1.8mm
2511



Plugger, 1.3mm-
1.8mm
2580



Burnisher,
1.75mm-2.25mm
2517



Burnisher & Plugger,
1.9mm-2.3mm
2512



Burnisher & Spatel,
1.5mm-2mm
2510



Composite Trimmer,
1.5mm-2.0mm
2557



Gingival Margin
Trimmer,
1.2mm, Mesial
2566



Gingival Margin
Trimmer, 1.2mm, Distal
2565



Plugger,
1.0mm-1.3mm
2579



Burnisher & Plugger,
Ballended,
1.5mm-2.5mm
2509



Packer-
Modeller IO
2575



Modeller Mh
2576



ErgoSteel

ErgoSteel Mini (S)

Length 100mm
Ø6.0mm

ErgoSteel Standard (SS)

Length 106mm
Ø8.0mm

ErgoSteel Plus (SP)

Length 106mm
Ø9.5mm

ErgoX

ErgoX

Length 100mm
Ø11.5mm

ErgoX Plus

Length 100mm
Ø13.6mm

Pluggers & Burnishers



**Modeller MH Mini
2577**



Burnisher



ErgoSteel

ErgoSteel Mini (S)

Length 100mm
Ø6.0mm

ErgoSteel Standard (SS)

Length 104mm
Ø8.0mm

ErgoSteel Plus (SP)

Length 106mm
Ø9.5mm

ErgoX

ErgoX

Length 100mm
Ø11.5mm

ErgoX Plus

Length 100mm
Ø13.6mm

Cement & Composite Spatulas



**Composite Spatula 6,
2.5mm-2.5mm**

2551



**Cement Spatula,
5mm-5mm
2568**



**Cement Spatula,
6mm-6mm
2569**



**Cement Spatula,
7mm-7mm
2570**



**Cement Spatel
2567**



**Mini Spatula, 2.0mm-
2.0mm
2550**



**Composite
Instrument Posterior,
2.2mm-2.5mm
2552**



**Universal Spatula,
1.5mm-2.0mm**

2571



Carvers



Nystrom I
2524

XP X SP SS S



Nystrom III
2525

XP X SP SS S



Heidemann, 2.0mm
2519

XP X SP SS S



Heidemann, 2.5mm
2520

XP X SP SS S



Carver 5
2523

XP X SP SS S



Heidemann, 3.0mm
2521

XP X SP SS S



Knife-Carver
2526

XP X SP SS S

Base Carriers



Base Carrier - Rigid
2506

XP X SP SS S



Base Carrier - Condenser Malleable Tip
2507

XP X SP SS S



Base Carrier - Condenser Rigid
2508

XP X SP SS S

Condensers



Amalgam Condenser,
1.0mm-1.5mm

2504



Amalgam Condenser,
2.0mm-2.5mm
2505



Filling
Contour
2564



Condenser, Drop-
Shaped,
1.3mm-1.8mm
2572



Form Condenser,
1.5mm-2.5mm
2562



Flow Instrument IO
2563



ErgoSteel

ErgoSteel Mini (S) Length 100mm Ø6.0mm

ErgoSteel Standard (SS) Length 104mm Ø8.0mm

ErgoSteel Plus (SP) Length 106mm Ø9.5mm

ErgoX

ErgoX Length 100mm Ø11.5mm

ErgoX Plus Length 100mm Ø13.6mm

Filling instruments



Filling Remover

2547



UP1 Plastic filling
2546



Filling
Instrument-C49
2501



Filling
Instrument-PFI 8A
2585



Filling Instrument-
PFI 6
2584



Filling Instrument-
PFI UP1
2583



Filling Instrument-
UP1
2582



Filling instruments



Filling Applicator

2531



XP X SP SS S



Filling Shaper

2534



XP X SP SS S




Filling Condensor

2535



XP X SP SS S



Special Filling Instrument,
2.0mm-2.5mm

2542



Plastic Instrument-Posterior

2548



Plastic Instrument-Anterior

2549



Cord Packers



Retraction Cord
Paker-UP1

2545



XP X SP SS S



Cord Packer, 2.0mm

2543



XP X SP SS S



Cord Packer, 2.6mm

2544



XP X SP SS S

ErgoSlip

Available in blue, gold, and black coatings with ErgoX, ErgoX Plus, and ErgoSteel Plus handles. Each set includes an ErgoTray-1 for secure handling.

Perfect For: Composite placement, restorative procedures, and general dentistry.

**Smoother
Restorations,
Every Time.**

Each ErgoSlip comes with an ErgoTray-1.



Burnisher &
Plugger, Ballended,
1.5mm-2.5mm

Plugger, 1.0mm-
1.3mm



2579K (Black)



2579G (Gold)



2579B (Blue)

Plugger, 1.3mm-
1.8mm



2580K (Black)



2580G (Gold)



2580B (Blue)

Plugger, 2.5mm-
3.8mm



2581K (Black)



2581G (Gold)



2581B (Blue)

Plugger, Cylindrical,
1.8mm-2.5mm



2578K (Black)



2578G (Gold)



2578B (Blue)

Burnisher



2513K (Black)



2513G (Gold)



2513B (Blue)

Burnisher & Plugger,
Ballended, 1.5mm-2.5mm



2509K (Black)




2509G (Gold)



2509B (Blue)

ErgoSlip

Smoother restorations,
Every Time.



Packer-Modeller IO



2575K (Black)



2575G (Gold)



2575B (Blue)

Modeller Mh



2576K (Black)



2576G (Gold)



2576B (Blue)

Modeller MH Mini



2577K (Black)



2577G (Gold)



2577B (Blue)

Knife-Carver



2526K (Black)



2526G (Gold)



2526B (Blue)

Universal Spatula,
1.5mm-2.0mm



2571K (Black)



2571G (Gold)



2571B (Blue)

Mini Spatula,
2.0mm-2.0mm



2550K (Black)



2550G (Gold)



2550B (Blue)

Condenser, Drop-
Shaped, 1.3mm-
1.8mm



2572K (Black)



2572G (Gold)



2572B (Blue)

Filling Instrument-
C49



2589K (Black)



2589G (Gold)



2589B (Blue)

Filling Instrument-
PFI 8A



2585K (Black)



2585G (Gold)



2585B (Blue)

Filling Instrument-
PFI 6



2584K (Black)



2584G (Gold)



2584B (Blue)

Filling Instrument-
PFI UP1



2583K (Black)



2583G (Gold)



2583B (Blue)

Plastic Instrument-
Anterior



2549K (Black)



2549G (Gold)




2549B (Blue)

Restoration

Non-stick composite instruments

ErgoSlip

Smoother restorations,
Every Time.



Plastic Instrument- Posterior



2548K (Black)



2548G (Gold)



2548B (Blue)

Filling Instrument- UP1



2582K (Black)



2582G (Gold)



2582B (Blue)

Filling Applicator



2531K (Black)



2531G (Gold)



2531B (Blue)

Filling Shaper



2534K (Black)



2534G (Gold)



2534B (Blue)

Filling Condensor



2535K (Black)



2535G (Gold)



2535B (Blue)

Filling Contour



2564K (Black)



2564G (Gold)



2564B (Blue)

Composite kits

All Round ErgoSlip (Non- Stick) Compo Set

9910

Available in



Instruments

Universal Spatula
2571B

Mini Spatula
2550B

Plastic Instrument, Anterior
2549B

Burnisher
2513B

Burnisher & Plugger 1.5 & 2.5mm
2509B


5 instruments + cassette

Composite kits

Super Restoration Kit

9912

Available in



Instruments

Mirror Handle with mouth Mirror
1401

Excavator 1.5mm
2802

Explorer 23 - Probe CP12
ball end
1311

Burnisher & Plugger 1.5
& 2.5mm
2509

Micro Sickle Scaler 01-02
2181

Mini Spatula
2550

Flow Instrument IO
2562


Tweezer Meriam, 15cm
3636-1A

8 instruments + cassette

Restoration Kits Compo Set 5

9911

Available in



Instruments

Form Condenser
2562

Mini Spatula
2550

Plastic Instrument Anterior
2549

Gracey 11-12
2563

Gracey 13-14
2182

5 instruments + cassette

Rubberdam instruments

ErgoDenta offers rubber dam clamps in three types:

- Gold-coated
- Black-coated
- Standard stainless steel



Gold-coated

Gold-coated rubber dam clamps offer superior corrosion resistance and a premium, non-reactive surface.



Standard steel

Standard stainless steel rubber dam clamps offer reliable durability and performance.



Black-coated

Black-coated rubber dam clamps provide reduced glare for improved visibility during procedures.

Clamps

Anteriors

00



2701-00
K2701-00
G2701-00

W00



2701-W00
K2701-W00
G2701-W00

6



2701-6
K2701-6
G2701-6

9



2701-9
K2701-9
G2701-9

9S



2701-9S
K2701-9S
G2701-9S

210



2701-210
K2701-210
G2701-210

212



2701-212
K2701-212
G2701-212

212SA



2701-212SA
K2701-212SA
G2701-212SA

106R



2701-106R
K2701-106R
G2701-106R

106L



2701-106L
K2701-106L
G2701-106L

B4



2701-B4
K2701-B4
G2701-B4

Lower Molars

3



2702-3
K2702-3
G2702-3

W3



2702-W3
K2702-W3
G2702-W3

7



2702-7
K2702-7
G2702-7

7A



2702-7A
K2702-7A
G2702-7A

W7



2702-W7
K2702-W7
G2702-W7

12A



2702-12A
K2702-12A
G2702-12A

W12A



2702-W12A
K2702-W12A
G2702-W12A

13A



2702-13A
K2702-13A
G2702-13A

26



2702-26
K2702-26
G2702-26

28




2702-28
K2702-28
G2702-28

202



2702-202
K2702-202
G2702-202

Pre molars




Using Rubber Dam Clamps Effectively

- **Select the Right Clamp:** Choose the correct size and shape for the tooth to ensure a secure fit.
- **Inspect and Sterilize:** Check for damage and ensure the clamp is clean before use.
- **Proper Placement:** Use forceps to position the clamp below the tooth's height of contour.
- **Check Stability:** Ensure the clamp is secure and won't dislodge during the procedure.
- **Ensure Comfort:** Adjust to avoid gingival impingement and maintain patient comfort.



Upper Molars




Molars / Universal




Clamps

Molars / Universal



Rubberdam instruments



Rubber Dam Kit

All-in-one kit with everything you need for quick and accurate rubber dam placement. The kit is available with standard, gold or black clamps

Instruments in the kit

Rubberdam clamps

2703-W8A	Universal: W8A
2705-W8	Universal: W8
2701-212	Anterior: #212
2701-9	Anterior molars: #9
2704-W2	Premolars: W2
2704-2	Premolars: #2
2704-2A	Premolars: 2A
2702-W7	Lower jaw molars: W7
2701-00	Lower jaw molars #00
2702-W3	Lower jaw molars: W3
2702-12A	Lower jaw molars : 12A
2705-30	Upper jaw molars: #30

Rubberdam instruments & accessories

2741	Ainsworth Rubber Dam Punch 17cm
2711	Brewer Rubber Dam Clamp Forcep 17cm
2761	Rubber Dam Clamp Organizing Board
2750	Rubber Dam Stand Small
2753	Rubber Dam Stand Large
5555	Sterilization cassette



9917
with standard clamps

9917G
with gold clamps

9917K
with black clamps

Crown instruments

Crown Scissors



Crown spreader



Crown remover pliers / forceps




Crown Remover incl. silicone tips
2953



Crown remover

Includes 3 tips, a wrench, and an accessory.
2951



Matrix retrainers

Nystorm



Tofflemire



Ivory



Periodontics



Gracey Curettes

ErgoDenta offers Gracey Curettes in Standard, Rigid, and Mini types, with ErgoX and ErgoSteel handles for comfort and durability.

Gracey Curettes



Gracey 1-2
2001

XP X SP SS S



Gracey 3-4
2002

XP X SP SS S



Gracey 5-6
2003

XP X SP SS S



Gracey 7-8
2004

XP X SP SS S



Gracey 9-10
2005

XP X SP SS S



Gracey 11-12
2006

XP X SP SS S



Gracey 13-14
2007

XP X SP SS S



Gracey 15-16
2008

XP X SP SS S



Gracey 17-18
2009

XP X SP SS S



Gracey 11-14
2010

XP X SP SS S



Gracey 12-13
2011

XP X SP SS S

Rigid Gracey Curettes



Rigid Gracey 11-12
2253



Rigid Gracey 13-14
2254



Gracey P3-P4
2014



ErgoSteel

ErgoSteel Mini (S)

Length 100mm
Ø6.0mm

ErgoSteel Standard (SS)

Length 104mm
Ø8.0mm

ErgoSteel Plus (SP)

Length 106mm
Ø9.5mm

ErgoX

ErgoX

Length 100mm
Ø11.5mm

ErgoX Plus

Length 100mm
Ø13.6mm

Mini Gracey Curettes



Mini Gracey 15-16
2106



Mini Gracey 7-8
2103



Mini Gracey 11-12
2104



Mini Gracey 13-14
2105



Mini Gracey 1-2
2101



Mini Gracey 17-18
2107



Sickle Scalers



Sickle Scaler H6-H7
2153



Sharp Jack Scaler
2154



Nebraska 128
2155



Sickle Scaler 204SD
2157



M23
2159



M23A
2161



Sickle Scaler 204S
2156



Crane-Kaplan
2162



Furcator Scaler
2168



Micro Sickie 01-02
2181



Mini Sickie 11-12
2182



Mini Excess Scaler
2184



ErgoSteel

ErgoSteel Mini (S) Length 100mm Ø6.0mm

ErgoSteel Standard (SS) Length 104mm Ø8.0mm

ErgoSteel Plus (SP) Length 106mm Ø9.5mm

ErgoX

ErgoX Length 100mm Ø11.5mm

ErgoX Plus Length 100mm Ø13.6mm

Universal Curettes



Columbia 2L-2R
2201



Columbia 4L-4R
2202



Columbia 13-14
2203



Barnhart 1-2
2204



Barnhart 5-6
2205



Langer 1-2
2206



Langer 3-4
2207



Langer 5-6
2208



Mini Langer 1-2
2209



Mini Langer 3-4
2210



Mini Langer 5-6
2211



McCall 13S-14S
2212



Mini McCall 13S-14S
2213



ErgoSteel

ErgoSteel Mini (S) Length 100mm Ø6.0mm

ErgoSteel Standard (SS) Length 104mm Ø8.0mm

ErgoSteel Plus (SP) Length 106mm Ø9.5mm

ErgoX

ErgoX Length 100mm Ø11.5mm

ErgoX Plus Length 100mm Ø13.6mm

Universal Curettes



**McCall 17-18
2214**



**McCall 17S-18S
2215**



More delicate and narrower
tip compared to 17-18.



**Duplex Curette
2217**



**Mini Duplex Curette
2218**



**Duplex Anterior
2219**



**Mini Duplex Anterior
2220**



**Goldman-Fox 3
2221**



ErgoSteel

ErgoSteel Mini (S)

Length 100mm
Ø6.0mm

ErgoSteel Standard (SS)

Length 106mm
Ø8.0mm

ErgoSteel Plus (SP)

Length 106mm
Ø9.5mm

ErgoX

ErgoX

Length 100mm
Ø11.5mm

ErgoX Plus

Length 100mm
Ø13.6mm

Hoe Scalers



**Hoe Scaler, Anterior
2281**



**Hoe Scaler, Lateral
2282**



**Hoe Scaler, Posterior
2283**



ErgoRazor

All ErgoRazors are available with ErgoX and ErgoX Plus handles. They come in black and gold coating or a silver finish without coating. The tips are made from 440C steel and undergo special cryogenic treatment in the UK for long-lasting durability. Each ErgoRazor comes with an ErgoTray-1 for efficient handling of the instruments.

Upto **3X** longer sharpness



Each ErgoRazor comes with an ErgoTray-1.



Gracey curettes

Gracey 1-2



2001RK (Black)

2001RG (Gold)

2001R

Gracey 3-4



2002RK (Black)

2002RG (Gold)

2002R

Gracey 5-6



2003RK (Black)

2003RG (Gold)

2003R

Gracey 7-8



2004RK (Black)

2004RG (Gold)

2004R

Gracey 9-10



2005RK (Black)

2005RG (Gold)

2005R

Gracey 11-12




2006RK (Black)

2006RG (Gold)

2006R

ErgoRazor

Upto **3X** longer sharpness



Gracey 13-14



2007RK (Black)

2007RG (Gold)

2007R

Gracey 15-16



2008RK (Black)

2008RG (Gold)

2008R

Gracey 11-14



2010RK (Black)

2010RG (Gold)

2010R

Gracey 12-13



2011RK (Black)

2011RG (Gold)

2011R

Mini Gracey 1-2



2101RK (Black)

2101RG (Gold)

2101R

Mini Gracey 7-8



2103RK (Black)

2103RG (Gold)

2103R

Mini Gracey 11-12



2104RK (Black)

2104RG (Gold)

2104R

Mini Gracey 13-14



2105RK (Black)

2105RG (Gold)

2105R

Sickle Scaler

Mini Gracey 15-16



2106RK (Black)

2106RG (Gold)

2106R

Sickle Scaler H6-H7



2153RK (Black)

2153RG (Gold)

2153R

Mini Sickle 11-12



2182RK (Black)

2182RG (Gold)

2182R


Micro Sickle 01-02



2181RK (Black)

2181RG (Gold)

2181R



Sickle Scalers

M23 Sickle Scaler

Ideal for efficiently removing supragingival calculus while maintaining clear visibility and precise control.

2159R
2159RK (Black)
2159RG (Gold)



XP X



Sharp Jack Scaler



XP X



2154RK (Black)

2154RG (Gold)

2154R

Nebraska 128



XP X



2155RK (Black)

2155RG (Gold)

2155R

Sickle Scaler 204S



XP X



2156RK (Black)

2156RG (Gold)

2156R

Sickle Scaler 204SD



XP X



2157RK (Black)

2157RG (Gold)

2157R

Universal curettes

M23A



XP X



2161RK (Black)

2161RG (Gold)

2161R

Columbia 2L-2R



XP X



2201RK (Black)

2201RG (Gold)

2201R

Columbia 4L-4R



XP X



2202RK (Black)

2202RG (Gold)

2202R

Columbia 13-14



XP X




2203RK (Black)

2203RG (Gold)

2203R

ErgoRazor

Upto **3X** longer sharpness



Barnhart 1-2



2204RK (Black)



2204RG (Gold)



2204R

Barnhart 5-6



2205RK (Black)



2205RG (Gold)



2205R

Langer 1-2



2206RK (Black)



2206RG (Gold)



2206R

McCall 13S-14S



2212RK (Black)



2212RG (Gold)



2212R

Mini McCall 13S-14S



2213RK (Black)



2213RG (Gold)



2213R

McCall 17-18



2214RK (Black)



2214RG (Gold)



2214R

McCall 17S-18S



2215RK (Black)



2215RG (Gold)



2215R

Duplex Curette



2217RK (Black)



2217RG (Gold)



2217R

Mini Duplex Curette



2218RK (Black)



2218RG (Gold)



2218R

Duplex Anterior



2219RK (Black)



2219RG (Gold)



2219R

Mini Duplex Anterior



2220RK (Black)



2220RG (Gold)



2220R

Goldman-Fox 3



2221RK (Black)



2221RG (Gold)




2221R

Perio Kits

Basic Hygienist Kit 9905

Available in



Instruments

Explorer 17-23
1107

Mirror Handle Special Blue
1401

Mini McCall 13s-14s
2213


Periodontal Probe CP12
1053

Mini Sickle Scaler
2182

5 instruments + cassette

All Round Perio Set 9903

Available in



Instruments

Mini Sickle Scaler
2182

Micro Sickle Scaler
2181

McCall 17-18 curette
2214

Gracey 1-2
2001

Gracey 9-10
2005

Gracey 11-12
2006

Gracey 13-14
2007

Gracey 15-16
2008


8 instruments + cassette

Basic Perio Kit 9906

Available in



X XP SP



Instruments

Mini Sickle Scaler
2182

Curette McCall 17-18
2214

Gracey 1-2
2001

Gracey 11-12
2006

Gracey 13-14
2007


5 instruments + cassette

Perio Diagnostic Kit 9907

Available in



X XP SP



Instruments

Micro Sickle Scaler
2162

Mirror Handle with Front
Surface Mirror No.4
1401

Explorer 8- Periodontal
Probe 2mm
1303

Gracey 11-12
2006

Gracey 13-14
2007


5 instruments + cassette

Tartar Remover Set 9904

Available in



X XP SP



Instruments

Explorer 17-23
1107

Sickle Scaler H6-H7
2153


Mini sickle scaler
2182

M23 Sickle Scaler
2159

Mirror Handle with Front
Surface Mirror No.4
1401

5 instruments + cassette

Extraction



Extraction Forceps

Upper Incisors, molars,
Premolars & Roots



Fig.1
Upper Incisors And Molars

S: 5005-1 | D:5025-1



Fig.17
Upper Molars, Right
S: 5005-17 | D: 5025-17



Fig.18
Upper Molars, Left
S: 5005-18 | D: 5025-18



Fig.30S
Upper Front
D: 5025-30S



Fig.29
Upper Front
S: 5005-29 | D: 5025-29



Fig.51A
Upper Roots
S: 5005-51A | D: 5025-51A



Fig.97
Upper Roots
S: 5005-97 | D: 5025-97



Fig.41
Upper Roots
S: 5005-41 | D: 5025-41 |



Fig.46
Upper Roots
D: 5025-46
S: 5005-46L (Very fine)
S: 5005-46XL (Very fine)



Fig.7
Upper Premolars
S: 5005-7 | D: 5025-7



Fig.67A
Upper Molars-For Upper Wisdom Teeth
S: 5005-67A | D: 5025-67A



Fig.OK
Universal-Upper Roots Extracting Forceps
S: 5005-OK



Diamond D



Standard S

Extraction Forceps

Upper Incisors, molars,
Premolars & Roots



Fig.98
Witzel (Upper roots)
D: 5025-98



Fig.144
Upper Roots
D: 5025-144



Diamond D



Standard S

Lower Incisors, molars,
Premolars & Roots



Fig.13.
**Lower Molars And
Premolars**
S: 5005-13 | D: 5025-13



Fig.22.
Lower Molars
S: 5005-22 | D:5025-22



Fig.23.
Lower Molars, Right
5005-23



Fig.24
Lower Molars, Left
5005-24



Fig.79
Lower Wisdom Teeths
S: 5005-79 | D: 5025-79



Fig.86A
Lower Molars
S: 5005-86A | D: 5025-86A



Fig.86B
Lower Molars
S: 5005-86B | D: 5025-86B



Fig.99
Hertz (Lower Molars)
D: 5025-99



Fig.74N
Lower Roots
S: 5005-74N | D: 5025-74N



Fig.UK
**Universal-Lower Roots
Extracting Forceps**
S: 5005-UK


Lower Incisors, molars, Premolars & Roots



Fig.87
Lower Molars
S: 5005-87



Fig.99
Hertz
D: 5025-99



Extraction forceps kit

Available in regular & diamond

Instruments in the kit	Diamond 5025	Regular 5005
Fig. 1. Upper Incisors And Molars	5025-1	5005-1
Fig. 13. Lower Molars And Premolars	5025-13	5005-13
Fig. 17. Upper Molars, Right	5025-17	5005-17
Fig. 18. Upper Molars, Left	5025-18	5005-18
Fig. 22. Lower Molars	5025-22	5005-22
Fig. 33A. Lower Incisors, Premolars & Roots	5025-33A	5005-33A
Fig. 51A. Upper Roots	5025-51A	5005-51A
Fig. 67A. Upper Molars-Upper Wisdom Teeth	5025-67A	5005-67A
Fig. 7. Upper Premolars	5025-7	5005-7
Fig. 79. Lower Wisdom Teeth	5025-79	5005-79

Diamond tip with a **serrated pattern** ensures a strong grip.



Diamond extraction
forceps set
5025



Standard extraction
forceps set
5005



10
Forceps
/box

Child Diamond Extraction Forceps



Fig.3
Upper Molar
S: 5010-3 | D: 5020-3



Fig.137
Upper Incisor
S: 5010-137 | D: 5020-137



Fig.139
Upper Pre Molar
S: 5010-139 | D: 5020-139



Fig.51S
Upper Roots
S: 5010-51S | D: 5020-51S



Fig.6
Lower Molar
S: 5010-6 | D: 5020-6



Fig.5
Lower Incisiver
S: 5010-5 | D: 5020-5



Fig.7
Lower Roots
S: 5010-7 | D: 5020-7

Child Extraction Forceps kit

Available in standard & diamond

Instruments in the kit	Diamond 5020	Regular 5010
Fig. 137. Upper Incisor	5020-137	5010-137
Fig. 139. Upper Premolar	5020-139	5010-139
Fig. 3. Upper Molar	5020-3	5010-3
Fig. 5. Lower Incisiver	5020-5	5010-5
Fig. 51S. Upper Roots	5020-51s	5010-51s
Fig. 6. Lower Molar	5020-6	5010-6
Fig. 7. Lower Roots	5020-7	5010-7



Diamond extraction forceps set 5020



Standard extraction forceps set 5010



Junior Extraction Forceps



Fig.22K
Molars, Lower Jaw
5015-22K



Fig.13K
Premolars, Lower Jaw
5015-13K



Fig.7K
**Upper Jaw-Molars
& Incisors**
5015-7K



Fig.38
**Incisors & Canins,
Lower Jaw**
5015-38



Fig.37
Incisors, Upper Jaw
5015-37



Diamond D



Standard S

Grip Master Diamond Forceps

ErgoDenta's GripMaster Extraction Forceps feature a diamond-coated beak surface, ensuring a secure hold with less grip strength required.



Upper Spade Root
D: 5040-2079



Upper Universal
D: 5040-2007



Upper Anterior
D: 5040-2076



Upper Molar
D: 5040-2021



Lower Universal
D: 5040-2002



Lower Anterior
D: 5040-1251



Lower Molar
D: 5040-2074



Lower Spade Root
D: 5040-145

Root tip elevators

ErgoDenta offers a variety of elevators, including Luxating Elevators, XTool, ErgoLite Elevators, and more, designed for comfort and long-lasting use.

Luxating Elevators



Straight



Name	Standard Handle	Classic Handle
Straight, 2.5mm	6001-1	6002-1
Straight, 3mm	6001-2	6002-2
Straight, 4mm	6001-3	6002-3
Straight, 5mm	6001-4	6002-4

Straight Gold



Name	Standard Handle	Classic Handle
Straight, 2.5mm	6001-1G	6002-1G
Straight, 3mm	6001-2G	6002-2G
Straight, 4mm	6001-3G	6002-3G
Straight, 5mm	6001-4G	6002-4G

Curved



Name	Standard Handle	Classic Handle
Curved, 2.5mm	6001-5	6002-5
Curved, 3mm	6001-6	6002-6
Curved, 4mm	6001-7	6002-7
Curved, 5mm	6001-8	6002-8

Curved Gold



Name	Standard Handle	Classic Handle
Curved, 2.5mm	6001-5G	6002-5G
Curved, 3mm	6001-6G	6002-6G
Curved, 4mm	6001-7G	6002-7G
Curved, 5mm	6001-8G	6002-8G

Luxating elevators kit

6002

Kit - 8 luxating elevators in classic handle



6002 G

Kit - 8 luxating Gold elevators in classic handle



Root tip elevators

Widman Knives



Straight
6090



Curved
6091

Bernard Elevators



Straight
6021-1



Angled
6021-2

Cryer



Right
6042-2



Left
6042-3

Appexo Elevators



Straight
2.5X9mm
6022-1



Straight
2.7mm
6022-2



Straight,
3mm
6022-3



Curved,
2.7mm
6022-5



Curved,
3mm
6022-6



Curved,
2.5mm
6022-4

Kopp Elevators



Edge Sharp,
3.6mm
6024-1



Left, 4.2mm
6024-5



Right, 4.2mm
6024-6



Round Sharp,
3.6mm
6024-2



Edge Sharp,
4.2mm
6024-3



Round Sharp,
4.2mm
6024-4

Root tip elevators

Lindo-Levian Elevators



Fig. 1
2.5mm
6046-1

Fig. 2,
3.5mm
6046-2

Fig. 3,
4.0mm
6046-3

Bein Root Elevators



Fig. 1,
2.3mm
6047-1

Fig. 2,
3.2mm
6047-2

Fig. 3,
4.1mm
6047-3

Flohr Elevators



Right
3.0mm
6045-2

Left
3.0mm
6045-3

Apical Root Elevators



Fig. 301, 3.0mm
6050-9

Fig. 302, 3.2mm
6050-10

Fig. 303, 3.3mm
6050-11

Hyllin Root Elevator



Fig. 1, 3.5mm
6053-1

Straight, 4mm
6053-2

Straight, Fig. 2,
4.5mm
6053-3

Heiderbrink Root Tip Elevators



Straight, Fig. 1,
2.5mm
6200-1

Right, Fig. 2,
2.5mm
6200-2

Left, Fig. 3,
2.5mm
6200-3

Root Tip Elevators




X-Tool elevators

Extract with Ease




X Tool-Luxating Elevators



6230

X Tool-Luxating Elevators Set

X Tool-Bernard Spade Elevators



X Tool-Bernard
Spade Tip
Elevators Set

6250



Root tip elevators

Berten Elevators



Straight
6023-1



Curved
6023-2



Angled
6023-3

Winter Elevators



Left Fig. 11L,
9.5mm
6080-1



Right Fig. 11R,
9.5mm
6080-2

ErgoVator

Meet ErgoVator. Luxating elevators with a color-coded design and ergonomic grip. Get ready for a smoother extraction experience—launching soon!



**Coming
soon**

Surgery



Scissors



Black Scissors

Experience sharp, accurate cutting with Super-cut scissors featuring a durable black coating.



- Minimal reflection
- Perfect for surgery
- Easy to clean
- Stays sharper longer
- Razor-like edge

Goldman-Fox



Straight
3007-4S
13 cm

Curved
3007-4C
13 cm

La-Garange



Double Curved
3008-4D
11.5 cm

IRIS



Straight
3002-4S
11.5cm

Curved
3002-4C
11.5 cm

Scissors



Standard

Reliable and versatile scissors for everyday cutting tasks.



TC

Tungsten carbide inserts ensure precise and long-lasting sharpness for delicate oral procedures.



Super-cut

High-performance scissors with micro-serrated blades for precise and clean cuts.

IRIS



Straight	Curved
3002-1S	3002-1C
11.5 cm	11.5 cm
Straight	Curved
3003-1S	3003-1C
15 cm	15 cm



Straight	Curved
3002-2S	3002-2C
11.5 cm	11.5 cm
3003-2S	3003-2C
15 cm	15 cm



Straight	Curved
3002-3S	3002-3C
11.5 cm	11.5 cm



Angular
3002-1A
11.5 cm

Scissors

La-Grange



Curved
3008-1D
11.5 cm



Curved
3008-2D
11.5 cm



Curved
3008-3D
11.5 cm

Goldman-Fox



Straight
3007-1S
13 cm

Curved
3007-1C
13 cm



Straight
3007-2S
13 cm

Curved
3007-2C
13 cm



Straight
3007-3S
13 cm

Curved
3007-3C
13 cm

Scissors

Metzenbaum Stump



Straight
3011-1S
14.5 cm



Curved
3011-1C
14.5 cm

Straight
3011-2S
14.5 cm

Curved
3011-2C
14.5 cm

Metzenbaum Sharp



Straight
3016-1S
14.5 cm

Curved
3016-1C
14.5 cm



Straight
3016-2S
14.5 cm

Curved
3016-2C
14.5 cm

Metzenbaum
Open



Straight
3013-2S
13.5 cm

Metzenbaum
Fino



Straight
3014-2S
14.5 cm

Curved
3014-2C
14.5 cm

Kelly Scissors



Straight
3026-1S
16 cm



Curved
3026-1C
16 cm



Straight
3026-2S
16 cm

Curved
3026-2C
16 cm

Scissors

Classic Suture Scissors



Straight 3009-1S 11.5 cm	Curved 3009-1C 11.5 cm
3010-1S 13 cm	3010-1C 13 cm

Kilner Scissor



Straight
3029-1S
12.5 cm

Joseph Scissor



Straight 3027-2S 14.5 cm	Curved 3027-2C 14.5 cm
---------------------------------------	-------------------------------------

Crown Scissor



Straight 3045-1S 13 cm	Curved 3045-1C 13 cm
-------------------------------------	-----------------------------------

Cooper Scissors



Straight
3052-1S
13 cm

3058-1S
16 cm
For Dressing

Spencer Scissors



Angled 3037-1A 11.5 cm	Straight 3039-1S 16 cm
-------------------------------------	-------------------------------------

Mayo Scissors



Straight
14.5 cm
3031-1S
3031-2S
18 cm
3032-1S
3032-2S

Curved
14.5 cm
3031-1C
3031-2C
18 cm
3032-1C
3032-2C

Micro Surgery Scissors

Perfect for small, delicate procedures

Noyes



Barraquer Scissors



Micro Surgery



- Straight 3107-1S 15 cm
- Curved 3107-1C 15 cm
- Offset 3107-1A 15 cm

- Straight 3121-2S 13 cm
- Straight 3045-2C 15 cm

- Straight 3101-2S 13 cm
- Curved 3101-2C 13 cm

Needle holders

Mayo Hegar



Straight
3212-1S
16 cm



Curved
3212-1C
16 cm



Straight
3212-2S
16 cm



Curved
3212-2C
16 cm



Straight
3226-1S
15 cm



Curved
3226-1C
18 cm

Fine
Swedish



Straight
3209-2S
15 cm

Crile
Ryde



Straight
3238-2S
15 cm

Crile-wood

TC needle holders offer superior grip and durability, ensuring precise and reliable suturing in surgical procedures.

Needle holders

Mathieu



Straight	Curved	Straight	Curved
3301-1S	3301-1C	3301-2S	3301-2C
14 cm	14 cm	14 cm	14 cm
3302-1S	3302-1C	3302-2S	3302-2C
17 cm	17 cm	17 cm	17 cm

Olsen



Straight
3312-2S
17 cm

Lichtenberg



Straight	Curved
3322-2S	3322-2C
17 cm	117 cm

Micro Needle holders

Castroviejo




Straight	Curved
3343-2S	3343-2C
16 cm	16 cm

Castroviejo
Micro



Straight	Curved
3331-2S	3332-2C
14 cm	18 cm



The Castroveijo Barraquer Needle Holder TC provides precise control and durability for detailed dental suturing.

Hemostatic Forceps

Halsted-Mosquito



Straight	Cruved
3407-1S	3407-1C
14.5cm	14.5 cm

Kocher
1x2



Straight	Cruved
3446-1SG	3446-1CG
13 cm	13 cm
3448-1SG	3448-1CG
16 cm	16 cm

Pean



Straight	Cruved
3411-1S	3411-1C
13.5 cm	13 cm
3413-1S	3413-1C
16.5 cm	16 cm

Crile



Straight	Cruved
3454-1S	3454-1C
14 cm	14 cm

Stieglitz



Curved 45°

Straight	Angled 45°
3456-1S	3456-1A
13 cm	13 cm

Removes splinters, broken files, metal posts, and small root fragments.

Allison



Straight
3432-1S
15 cm

Adson



Straight	Curved
3461-1S	3461-1C
14 cm	14 cm
3462-1S	3462-1C
18 cm	18 cm

Kelly



Straight	Angled
3441-1S	3441-1C
14 cm	14 cm

Hemostatic Forceps

Micro-Mosquito



Straight
3401-1S
10 cm
3402-1S
12 cm

Curved
3401-1C
10 cm
3402-1C
12 cm

Hemostatic Forceps



3437-1S
12 cm

Tohoku



3478
10 cm

Backhaus



3470 9 cm
3471 11 cm

Tweezers

College



Serrated
3627-1A
15 cm

Short Serrated, 13mm
3628-1A
15 cm



Serrated, with Lock
3630-1A
15 cm



Serrated
3636-1A
16 cm



Serrated
3636-2A
16 cm

Meriam

Surgical
1x2



Kocher
3633-1G
14.5 cm

Tweezers

Adson



3606-2
12 cm



3606-1
12 cm

1x2



Kocher
3606-1G
12 cm
Kocher
3608-1G
15 cm

offset



Kocher
3608-1AG
15 cm

Tweezers by ErgoDenta



Features serrated or smooth tips for precise grip.

Perry



3621-1C
13 cm

Flagg



Half Serrated
3616-1
16 cm

Anatomical



3646-1
13 cm
3647-1
14.5 cm
3648-1
16 cm



3646-2
13 cm

Miller



Serrated
3622-1
16 cm



Serrated
3622-1C
16 cm

Tweezers

Semken



3657-1
15 cm



3656-1
12 cm



3657-1G
15 cm



3657-1AG
15 cm



3656-1C
12 cm



3657-1C
15 cm

Endo Tweezers

DeBakey

Taylor

Crane Kaplan



Cross hold

3674-1A
15cm



With lock

3676-1A
15cm



3675-1A
15cm



Offset
3667-1
16 cm



3667-1A
16 cm



3651-1
15 cm



1x2
3651-1G
15 cm



3662-1AG
15 cm
Left



3661-1AG
15 cm
Right

Tweezers

Micro surgical tweezers

Gerald

Micro surgical plier



3801-1G
15 cm



Anatomical
3802-1S
18 cm



Curved
3802-1CG
18 cm



1x2
3671-1G
17 cm



Curved
3671-1CG
17 cm




3724-1C
15 cm

Lip, cheek & mouth retractor

Sean-Miller

Obwegeser
Internal

Obwegeser
External



4007
Blunt

4006
Sharp



Small
32x10mm
4053
Medium
42x10mm
4054
Large
55x12mm
4055

23cm



Small
14x17mm
4036
Medium
20x22mm
4037
Large
26x30mm
4038
21 cm

Lip, cheek & mouth retractor

Wassmund



Langenbeck



Middeldorpf



Bowdlerhenry's



2x16cm
15.5cm
4127
Cawood-
Minnesota



Fig.1, 26 x 10 mm / 30 x 13 mm
Fig.2, 30 x 10 mm / 34 x 13 mm

4101
Farabeuf
Set of 2 pieces
12cm



Lip, cheek & mouth retractor

2cm x 14cm

4106

Minnesota



4041

Universal
cheek
retractor

14cm



4042

Universal
cheek
retractor

14cm



Syringes

Self Aspirating Syringe

Cartridge Syringe



18ml

3929



22ml

3930



18ml

3931



18ml

3901



22ml

3902



22ml

3910



22ml

3907

Pistol model
syringes

Pen model
syringes

Syringe
holder

0.20 ml
per click
1.8ml

3961

Spare-parts
available



1.8ml

3968

Consists of three
detachable
components

0.07 ml/ click



1.8ml

3969

The offset/curved
model syringe is
versatile for use in
both the upper and
lower jaw.


0.07 ml/ click





3999

Gutta-Percha Remover



Hook, Remove & Move On!

The ultra-thin 0.4 mm tip preserves canal walls. The double-ended design offers versatile hooks and easy retrieval from any angle, leaving no residue.

7001

Gutta Percha Remover, Up-Down, 0.4mm Hook


7002

Gutta Percha Remover, Left-Right, 0.4mm Hook

Endo instruments


Endo Plugger

Used for condensing gutta-percha, providing a dense and secure canal filling.



Endo Spreader

Designed to laterally compact gutta-percha, creating space for additional filling material for a superior seal



ErgoSteel



ErgoX



Endodontic Essentials

Browse 600+ endodontic products, including the Endo Black series, instruments, accessories, kits, and more in our ErgoEndo catalog. Available soon!

Endo instruments

Endo Spreader & Pluggers

5 | 10 | 15 | 20



Endo Spreader & Pluggers
Ø0.4mm & Ø0.4mm

7056



Endo Spreader & Pluggers
Ø0.4mm & Ø0.6mm

7057



Endo Spreader & Pluggers
Ø0.6mm & Ø0.8mm

7058



Endo Spreader & Pluggers
Ø0.6mm & Ø1.0mm

7059



ErgoSteel


ErgoSteel Mini (S)

Length 100mm
Ø6.0mm



ErgoSteel Standard (SS)

Length 104mm
Ø8.0mm



ErgoSteel Plus (SP)

Length 106mm
Ø9.5mm



ErgoX

ErgoX

Length 100mm
Ø11.5mm



ErgoX Plus

Length 100mm
Ø13.6mm



Endo Explorers & Pluggers



Endo Explorer

7050



Root Canal
Explorer

7051



Endo Plugger,
0.5mm-1.0mm

7055



Endodontics Forceps



Endodontics
Forceps

7411



Endodontics File Holder
7415



Orthodontics



Separating Plier for an Easy Placement of Separating Elastics, 15cm
2730



Schwarz, Fig. 1, max Ø 0.4mm, 13 CM
4266



Adams Plier, max Ø 0.4mm, Lap joint, 13cm
4310



Universal Plier, TC
4382



Angle Plier, max Ø 0.7mm, Box Joint, 12.5cm
4312



Angle Plier, max Ø 0.7mm, Lap Joint, 12.5cm
4313



Wire Cutter, max Ø 1.0mm, 15cm
4320



Wire Cutter, max Ø 0.6mm, 12cm, TC
4321



Langback Plier, max Ø 0.8mm, 12cm
4344



Weingart Plier for Dental Technic, 15cm, TC
4273



Wire Cutter, max Ø 0.6mm, 12cm
4322



Micro ligature cutter
4379

Surgery

Orthodontics



Tweed Loop
Plier, 13cm
4294



Band Pusher
4372



Ligature Cutter
4275



Measuring
ruler, 0-60mm
4373



Measuring ruler,
0-80mm
4374



Cinch Back
4375



Weingart Plier Standard, for
Dental Technic, 15cm
4272



Boon's Gauge
With Lid
4377



Band seater,
Blue, 13cm
4390B



De la rosa band
contouring
4378



Pin & Ligature
Wirecutter, max Ø
0.5mm, 13cm, TC
4301



Pin & Ligature
Wirecutter, max Ø
0.7mm, 13cm
4302



Pin & Ligature
Wirecutter, max Ø
0.7mm, 13cm, TC
4303

Surgery

Orthodontics



Bracket
Tweezer, 11cm
3694-1A



Aderer Plier, Short.
max Ø 0.6mm,
11.5cm
4230



Aderer Plier,
Long, max Ø
0.6mm, 13cm
4231



HOW Plier,
Straight, 14cm
4232



Distal End Wirecutter
Max. Ø 0.55mm-Ø
0.64-12.5 CM
4227



20mm

Flat Nose Pliers,
Short, 15cm
4252



25mm

Flat Nose Pliers,
Long, 15cm
4253



Buccal Tube Holder
4376



Tweed, max Ø
0.5mm, Short, 12.5cm
4296



Hard wire cutter, 10°
Angle, 14cm, TC
4290



Ribbon-Tweed
arch, Short
4292



Ribbon-weed
arch, Long
4293



Three beak Plier
4223



Tweed Wire
Bending Plier, 13cm
4295

Surgery

Orthodontics



Adhesive Removal Plier
4201

Changeable Part for
Adhesive Removal
Plier, Blade
4201E 2

Changeable Part for
Adhesive Removal
Plier, Teflon
4201E 1



**Ortho Arch Wire
Tucker**
7072



Ligature 17C
7073



Ligature U15
7074



**Ortho Ligature
Tucker**
7071



**Bracket
Removal Plier**
4280



**Bracket Removing
Plier, 12.5cm**
4281



**Posterior Band Removing
Plier, Long, With
Teflon/Plastic Tip, 15cm, TC**

4206

Changeable Part for
Posterior Band
Removing Plier, Plastic
Pad

4206E1



**Mathieu Ortho
Plier**
4371



Mathieu Ortho plier
4370



**Distal End Cutter,
12cm, TC**
4225



**Bracket Positioning
Tweezer, 14cm**
3692-1

Surgery

Implantalogy

These impression trays with retentions rim are made from top-quality material, providing a firm grip that won't bend and ensuring consistent, accurate impressions.

Perforated Impression Tray

Upper Jaw



4521-1 XS



4521-2 S



4521-3 M



4521-4 L



4521-5 XL

Lower Jaw



4501-1 XS



4501-2 S



4501-3 M



4501-4 L



4501-5 XL

Non - Perforated Impression Tray

Upper Jaw



4571-1 XS



4571-2 S



4571-3 M



4571-4 L



4571-5 XL

Lower Jaw



4551-1 XS



4551-2 S



4551-3 M



4551-4 L



4551-5 XL

Half impression trays

Perforated



Non-Perforated



Impression trays Kits

Perforated set



4521

5 upper & 5 lower perforated impression trays set

Non-perforated set



4521

5 upper & 5 lower non perforated impression trays set

Surgery

Implantology

Sinus-lift instruments



All Sinus-lift instruments comes in ErgoSteel Plus handle

Curettes #1
7201

4mm-4mm



Curettes #2
7202

3mm-4mm



Kramer-Nevins #1
7206

4mm-7mm



Kramer-Nevins #2
7207

4mm-7mm



Sinus Lift Instrument
7211

4mm-4mm



Sinus Lift Instrument
7212

3mm-3mm



Sinus Lift Instrument
7213

4mm-4mm



Sinus Lift Instrument
7216

4mm-4mm



Sinus Lift Instrument
7217

4mm-4mm



Sinus Lift Instrument
7218

4mm-4mm



Hammers



18cm

Spare jaws for mead mallet

4401
Teflon hammer
Ø26mm



17cm

4402
Solid head hammer
Ø20mm



19cm

4403
Solid head hammer
Ø20mm

Implantology kits

Osteotome

Osteotome straight

Set of 5 including hammer & cassette

7270

Osteotome size

- 2.6 mm
- 3.1 mm
- 3.8 mm
- 4.5 mm
- 5.0 mm



Osteotome offset

Set of 5 including hammer & cassette

7271

Osteotome size

- 2.6 mm
- 3.1 mm
- 3.8 mm
- 4.5 mm
- 5.0 mm



Condenser

The condenser kit ensures each tip fits the previous cavity, allowing safe depth measurement and matching standard implant lengths.

Condenser straight

Set of 7 including hammer & cassette

7280

A	B
2mm	2.7mm
2.5mm	3.2mm
3mm	3.7mm
3.5mm	4.2mm
4mm	4.7mm
4.5mm	5.2mm
5mm	5.7mm



Condenser offset

Set of 7 including hammer & cassette

7281

A	B
2mm	2.7mm
2.5mm	3.2mm
3mm	3.7mm
3.5mm	4.2mm
4mm	4.7mm
4.5mm	5.2mm
5mm	5.7mm



Bone Spreader

Set of 6 including hammer & cassette

7282

A	B
16mm	16mm
16mm	22mm
22mm	22mm



Bone Splitting

Set of 8 including cassette

7250

Straight	Curved
3mm	3mm
4mm	4mm
6mm	6mm
8mm	8mm



Bone curettes

Lucas



Lucas 84
 2.0mm
 7103



X XP SP SS S



Lucas 85
 2.5mm
 7104



X XP SP SS S



Lucas 86
 3.0mm
 7105



X XP SP SS S



Lucas 87
 3.5mm
 7106



X XP SP SS S



Lucas 88
 4.0mm
 7107



X XP SP SS S

Hemingway



Hemingway
 2mm
 7110



X XP SP



Hemingway-2
 2.8mm
 7112



X XP SP



Hemingway-3
 3.5mm
 7113



X XP SP

Hemingway Straight



7114
 2.8 - 3.5mm

ErgoSteel

ErgoSteel Mini (S) Length 100mm Ø6.0mm
 ErgoSteel Standard (SS) Length 104mm Ø8.0mm
 ErgoSteel Plus (SP) Length 106mm Ø9.5mm

ErgoX


ErgoX Length 100mm Ø11.5mm
 ErgoX Plus Length 100mm Ø13.6mm

Bone curettes

Miller



Periosteal elevators




Surgery instruments


Bone files



Periotome



Bone rongeurs



Pro Surgery kit 9916

A complete set of precision surgical instruments, including micro tweezers, periosteal elevators, scissors, retractors, forceps, and essential accessories, all organized in a large cassette for efficient procedures.



7102 SS
Miller 21



3801-1G
Micro Surgery
Tweezer, 1x2,
15cm



2547 SS
Filling Remover



3007-4C
Goldman Fox
Scissor TC curved,
Super Cut, Black



1310 SP
Expros 23-Probe
CP 12



3209-2S
Fine Swedish
Needle Holder TC



4014
Langenbeck
Retractor, 42 x11
mm, 22cm



1450 SS
Mirror Handle,
8mm



7156 SP
Periosteal
Elevator, Molt



7158B SP
Glickmann,
titanium coated,
periosteal elevator



3437-1S
Hemostatic Forcep



5791
Scalpel Handle



3627-1A
College Tweezer,
15cm



5594
Bowl 50x60x15mm



3008-2D
La garange
Suture scissor
TC



3657-1
Semken Straight



1310 SP
Expros 23-Probe
CP 12



3999
Syringe Holder



3929
Self Aspirating
Syringe, Plastic &
Metal



5603
Bowl / Pots Ø 70
mm (190ml)



3801-1G
Large Cassette
with removeable
top for 20 pcs

All Round Surgery Kit 9913

A versatile set for various surgical procedures, including tweezers, retractors, scissors, probes, and essential instruments, all organized in a sterilization tray.




Micro Surgery Kit 9914

A specialized set for delicate procedures, featuring micro tweezers, periosteal elevators, precision scissors, a needle holder, and essential tools, all neatly organized in a sterilization tray.




Diverse



Sterilization Cassettes with detachable lid

Stainless steel cassettes are rust-free. The detachable-lid design allows for better instrument visibility, as well as saving space within the cassette, easy handling, and open/close functionality. Fully autoclavable and capable of disinfection as well as sterilization in a single cassette. Our trays are available in five different colors for easy identification: red, blue, green, orange and yellow.



Small tray for 5 instruments

5511

H:35mm X
W:145mm X
L:185mm

5511B
5511G
5511R
5511O
5511Y



Small tray for 10 instruments

5512

H:35mm X
W:145mm X
L:185mm

5512B
5512G
5512R
5512O
5512Y



Small tray for 20 instruments


5513

H:35mm X
W:145mm X
L:185mm

5513B
5513G
5513R
5513O
5513Y

Standard Sterilizations Cassettes

Simplify and reduce instruments processing time
Disinfect & sterilization of instrument in one tray



5531B
5531G
5531R
5531O
5531Y

Sterilization Cassettes for 5

5531

H:35mm X
W:145mm X
L:185mm



5532B
5532G
5532R
5532O
5532Y

Sterilization Cassettes for 7

5532

H:35mm X
W:145mm X
L:185mm



5533B
5533G
5533R
5533O
5533Y

Sterilization Cassettes for 10

5533

H:35mm X
W:145mm X
L:185mm



5534B
5534G
5534R
5534O
5534Y

Sterilization Cassettes for 20

5534


H:35mm X
W:145mm X
L:185mm



ErgoTrays

Daily Essentials for your practice

Key Features



ErgoTray 5



The ErgoTray 5 is ideal for sterilizing up to five dental instruments. It has a secure lock and an ergonomic lid that prevent tip injuries and ensure thorough cleaning.
180 x 86 x 26 mm



Orange 5505O



Pink 5505P



Red 5505R



Yellow 5505Y



Blue 5505B



Black 5505K




Grey 5505G



Purple 5505E

ErgoTrays 8



The ErgoTray 8 is ideal for sterilizing up to eight dental instruments. It has a secure lock and an ergonomic lid that prevent tip injuries and ensure thorough cleaning.
180 x 136 x 26 mm



Orange 5508O



Pink 5508P



Red 5508R



Purple 5508E



Yellow 5508Y



Blue 5508B



Black 5508K



Grey 5508G

ErgoTray 1

The ErgoTray single instrument holder is ideal for organizing and managing your ErgoX and X Plus instruments. This tray is included with all ErgoRazor and ErgoSlip instruments, and it is both fully autoclavable and reusable.



Grey 5501I



Black 5501K

ErgoVator Tray

The ErgoVator tray is the perfect solution for managing your ErgoVator elevators. Each ErgoVator elevator is supplied with an ErgoVator tray included.



Grey 5502I



Black 5502K

Solid & Perforated trays



Solid Tray
5582
180x280mm



Perforated Tray
5583
180x280mm



Perforated Base Tray
5584
180x280mm



Instrument Tray
5616
250mm x 150mm
x 40mm



Shallow (Mayo) Tray
5617
250mm x 150mm
x 40mm



Scaler Tray
5618
200mm x
100mm x 20mm



Perforated Tray
5620
200mm x 300mm



Catheter Tray
5621
220mm x 120mm x 50mm

**Storage Tubes/
Instrument Holder Ø**
5585
50x100mm



**Storage Tubes/
Instrument Holder Ø**
5587
50x170mm

**Storage Tubes/
Instrument Holder Ø**
5588
75x190mm

**Storage Tubes/
Instrument Holder Ø**
5589
90x210mm



Holders, boxes & bowls



**Instrument box
with Lid**
5590
100x200mm



Kidney Shape
5592
200mm



Bowl
5594
50mmx60mmx15mm

Bowl / Pots Ø 55 mm
5601
(65ml)

Bowl / Pots Ø 65 mm
5602
(130ml)

Bowl / Pots Ø 70 mm
5603
(190ml)



Round Box With Lid
5623
120mm x 120mm



Glass (18cl)
5624
70mmx80mm



**Instrument
Tray With Lid**
5622
220mm x 120mm x 50mm

Bowl / Pots Ø 70 mm
5604
(90ml)

Bowl / Pots Ø 90 mm
5605
(180ml)

Bowl / Pots Ø 110 mm
5606
(360ml)



Scalpel Handles



No:3 180°
5663



No:3
5664



No:3
5665

Accessories



Rulers Scalpel 10mm,
5661
12 pcs per pack



Rulers Scalpel 15mm,
5662
12 pcs per pack



Napkin Holder
5666
approx 44cm

Calipers



Sliding Calliper
5673
0-125mm



Iwanson
5674
10mm



Callipers Beerendonk
5667
130mm



Zuricher
5668
130mm



Boley
5669
145mm

Bur Holders

These bur holders can hold 4, 8, or 12 burs and are suitable for any FG or latch type. The colored silicone holders are replaceable and come in 4 colors, helping to keep things organized easily.



For 4 Bur Holder
7381B



For 8 Bur Holder
7382R



For 12 Bur Holder
7383G

Endo accessories

Fully autoclavable endo accessories available in various colors and sizes. View more on our website



Bur Disinfection box resistant to high temperature and Pressure, Blue
7458B



Translucent Bur disinfection box Square shape, 36 Holes, Blue
7454B



Gutta Percha Points Shaping Tip Cutter, Brown
7460N



Root Canal File Frame, 16 Holes, Blue
7459B



Endodontic block with memory holder, 5 holes, Orange
7462O



Gutta Percha Points Measuring Ruler, Black
7461K



Root Canal File Frame, 6 holes, Grey
7463I




Root Canal File Frame, 8 holes, Sky Blue
7453S




Endo File Cleaner Round Shape, Sky Blue
7452S



Translucent Bur disinfection box Rectangle shape, 72 Holes, Blue
7455B



Opening Disinfection Frame Resistant to high temperature and pressure, 25 Holes, Orange
7457O



Equipment of high temperature and high pressure disinfection box, 4 Holes, Green
7456G



Wrist Endo Measuring Cleaner, Sky Blue
7451S



72 holes endo box for diamond bur, Purple
7464M



Root canal measuring plate, Grey
7460G

File holder



Light Blue
7421I

Gray
7421B



Color-code rings

Organize and identify your instruments quickly and efficiently with our fully autoclavable, heat-resistant color-code rings. Available in two sizes—small for hand instruments and large for extraction forceps or instruments with bigger diameters—they ensure hassle-free color-coding and a secure fit.

Small:
50 pieces per box per color

Large:
25 pieces per box per color

- Fully autoclavable and heat-resistant
- Available in 8 vibrant colors

Color	Large	Small
Blue	5641B	5651B
Green	5642G	5652G
Red	5643R	5653R
Yellow	5644Y	5654Y
Black	5645K	5655K
White	5646W	5656W
Purple	5647P	5657P
Light blue	5648S	5658S

ErgoTape

3x longer than regular tapes

- 10-meter roll length
- Ø5mm thickness
- Available in 8 vibrant colors
- Fully autoclavable
- Built-in cutter for easy, fuss-free use



Color	Small
Blue	5691B
Light Yellow	5691L
Red	5691R
Dark Yellow	5691Y
Black	5691K
White	5691W
Purple	5691M
Light Pink	5691P
Assorted pack	5691

3x longer than regular tapes

Built in cutter for easy tape cutting

8 colors

Fully autoclavable

ErgoLite

Ergonomic, Durable & Cost-
Effective Instruments



ErgoLite X

ErgoLite X instruments are lightweight (11.5 mm, 14 g) and ergonomic in order to reduce hand strain. Made from anti-rust stainless steel, they are fully autoclavable and feature distinct colors and clear markings for easy identification and everyday use.





Diagnostics Explorers



Explorer + Probes

Explorer & Probe in One – Efficient & Convenient Diagnostics.



Endo instruments



Restoration Excavators

Action Excavator



Round Excavator



Filling instruments



ErgoLite X

ErgoLite X instruments are lightweight (11.5 mm, 14 g) and ergonomic in order to reduce hand strain. Made from anti-rust stainless steel, they are fully autoclavable and feature distinct colors and clear markings for easy identification and everyday use.



Restoration Filling instruments



Nystrom III
9120



Cord Packer 2.0 mm
9119



Burnisher & Plugger 1.5-2.5mm
9118



Burnisher 1.75-2.25 mm
9127



Form Condenser
9126



Plastic Instrument (Anterior)
9125



Composite Spatula 6
9124



Nystrom I
9123

Periodontics Mini Gracey



Mini Gracey 13-14
9129



Mini Gracey 11-12
9128

Periodontics Gracey curettes



Gracey 13-14
9132



Gracey 11-12
9131



Gracey 7-8
9133



Gracey 1-2
9130

Sickle Scalers



204SD
9140



204S
9139



H6-H7
9138



M23
9137



Crane-Kaplan
9136



Micro Sickle Scaler
9135



Mini Sickle Scaler
9134

ErgoLite S

ErgoLite Steel offers two handle options to suit your preference—solid steel for durability or a lightweight, ergonomic hollow handle for enhanced comfort. Both feature a classic textured design for a secure grip, making them perfect for everyday use.



Diagnostics

Explorers


Single-end

Double-end




Probes

Single-end



Explorer + Probes




Restoration

Excavators



Filling instruments



ErgoLite S

ErgoLite Steel offers two handle options to suit your preference—solid steel for durability or a lightweight, ergonomic hollow handle for enhanced comfort. Both feature a classic textured design for a secure grip, making them perfect for everyday use.



Restoration

Filling instruments




Periodontics

Gracey Curettes



Sickle Scalers



Bone curettes

Lucas



Hemmingway



Periosteal elevators


Molt



ErgoLite Elevators

Now Available: Gold Titanium-Coated Luxating Elevators!

Along with our regular luxating elevators, we now offer gold titanium-coated luxating elevators for enhanced durability, reduced wear, and longer-lasting sharpness for smoother extractions.



Comes in the tray included

Straight




Name	Regular	Gold
Straight, 1mm	9204	9204G
Straight, 2mm	9205	9205G
Straight, 3mm	9206	9206G
Straight, 5mm	9207	9207G

Curved



Name	Regular	Gold
Curved, 3mm	9201	9201G
Curved, 5mm	9203	9203G

Duel edge



Name	Regular	Gold
Duel edge, 3mm-1.5mm	9208	9208G
Duel edge, 5mm-3mm	9209	9209G

Contra/Inner Curved



Name	Regular	Gold
Inner Curved, 3mm	9202	9202G
Contra Curved, 3mm	9210	9210G



Thank You for Browsing Our Catalog

Explore More on Our Website

Maintenance of Dental Instruments

Learn how to properly care for your instruments to ensure long-lasting performance.

Warranty Information

Find out more about our warranty coverage and how we support our products.

Orders & Returns

Get all the details on how to place orders and handle returns quickly and easily.

Become a Distributor or Reseller

Interested in partnering with us? Visit our page to learn how to become an official distributor or reseller.

Private Labeling

Get everything you need—from design to final product—through our one-stop private label solution. No intermediaries, just higher ROI and total brand control.

ErgoDenta ApS
Hørskættø 14, 2630 Taastrup | Denmark
Tel. +45 5524 1100
www.ergodenta.dk | info@ergodenta.dk